

New And Efficient Method For Scaling Up Human Mesenchymal Stem Cells

¹Sanyal, Suparna, ¹Slater, Katie, ¹Abraham, Elizabeth, ¹Flaherty, Paula, ¹Medini, Tony, ¹Spargo, Cathy, ¹Reyes, Jamie, ¹Wongden, Eileen, and ¹Hastings, Abel
¹BD Biosciences, Ten Oaks Park, Bedford, MA USA
¹BD Technologies, 21 Davis Dr., Research Triangle Park, NC USA



Abstract

The use of single layer T-flasks can prove to be time consuming, cumbersome and labor intensive for cell culture applications requiring large number of cells. One solution needs to be exercised while scaling-up by using multi-layer vessels as differences in culture environment between layers of multi-layered vessels can affect yield, homogeneity and potential functional performance of cell populations. In the current study, we tested whether scaling-up human mesenchymal stem cells (hMSCs) in BD Falcon® Multi-Flasks affects their growth or characteristics compared to cells grown in control T-175 flasks. BD Falcon Multi-Flasks have been designed with an equalization/mixing port that enables uniform distribution of cells and media on each layer of the culture vessel thereby providing a homogeneous culture environment for cells throughout. hMSCs grown in media with and without serum, displayed no difference in morphology, attachment, proliferation, viability or following culture in BD Falcon Multi-Flasks compared to cells cultured in control T-175 flasks. To verify impact of long-term culture, hMSCs were thawed directly in both vessel types and cumulative population doublings was determined for five consecutive passages. hMSC expansion potential was found to be equivalent per unit area in both vessel types. The ability of cell passage kits to maintain multipotency and multipotency was also compared between cells cultured in single and multi-layered vessels. In addition to the hMSCs, diverse established and primary cell lines have also been cultured in the BD Falcon Multi-Flasks with equivalent performance to T-175 flasks. This study demonstrates a new and efficient way to scale-up stem cells without the need for re-optimizing existing cell culture conditions or compromising the quality, performance and homogeneity of cell populations.


Materials and Methods

Passaged Human Mesenchymal Stem Cells (Lonza) were cultured as per vendor's guidelines. Adipogenic and osteogenic differentiation kits were also purchased from Lonza. Oil-Red-O stain (Miles) was used to stain adipocytes and Vector Blue (Vector Laboratories, Inc.) was used to stain osteocytes positive for alkaline phosphatase activity. CellTrace™ CFSE fluorescent dye for proliferation assays was purchased from Invitrogen. hMSCs were cultured in either serum-containing growth medium (hMSC-GM/ Lonza) or in a serum-free medium (BD Falcon™ hMSC SF CCE). For serum-free culture, all vessels were coated with appropriate attachment factors. Cell growth and characterization studies using hMSC-GM and BD Falcon™ hMSC SF CCE media was conducted in independent laboratories to evaluate the performance of BD Falcon Multi-Flasks compared to control T-175 flasks.

Conclusions


- hMSC growth in media with and without serum was equivalent in BD Falcon Multi-Flasks and control T-175 flasks.
- hMSC morphology, attachment and proliferation rates were equivalent following culture in BD Falcon Multi-Flasks vs. T-175.
- hMSC cell expansion was higher and more rapid in BD Falcon™ hMSC SF CCE media compared to serum containing media.
- Immunophenotype of late passage hMSCs were identical between cells grown in single and multi-layer vessels.
- The extent of differentiation observed with late passage hMSCs following culture in BD Falcon Multi-Flasks was comparable to cells grown in control T-175 flasks.
- BD Falcon Multi-Flasks provide a uniform culture environment to efficiently scale up stem cells as well as other diverse cell types from T-flasks without the need for lengthy protocol re-optimization.

1 Equivalent Cell Growth and Cumulative Population Doublings of hMSCs in BD Falcon Multi-Flasks and Control T-175 Flasks




Cell yield per unit area and cumulative population doublings were equivalent for hMSCs grown in BD Falcon Multi-Flasks and T-175 flasks in media with and without serum. Four to five thousand cells/cm² were seeded in control T-175 and 3-layer BD Falcon Multi-Flasks in triplicate vessels and cultured in media with and without serum until approximately 80% confluency. Cells were harvested, counted and cell yield per unit area was determined. Cell yield per unit area for hMSCs after passage 1 cell growth. Cumulative population doublings of hMSCs from each vessel type were calculated for five consecutive passages (right panel).

2 Equivalent Cell Morphology, Attachment and Proliferation from BD Falcon Multi-Flasks and Control T-175 Flasks




A) Cell attachment and growth profile of hMSCs were equivalent for cells cultured in BD Falcon Multi-Flasks and control T-175 flasks. Cell attachment, spreading and growth was measured using an impedance-based assay (CELLGrowth™ RTCA MP System, Roche) following culture of hMSCs in hMSC-GM for 5 days in each vessel type. Cells were plated in serum-free media in 96 well 96-plates and interactions between cells and electrodes measured as changes in impedance, were reported as arbitrary cell index (CI) values. CI values are positively correlated with increased cell attachment and growth. B) Morphology of hMSCs cultured in BD Falcon Multi-Flasks were indistinguishable from those cultured in control T-175 flasks. C) Proliferation profile of hMSCs were very similar for cells cultured in both, single and multi-layered vessels. hMSCs were stained with 2.5 mM CellTrace™ CFSE fluorescent dye for proliferation (CFSE) for 19 minutes, washed and seeded in control T-175 and BD Falcon Multi-Flasks in hMSC-GM. Five days after culture, cells were harvested and analyzed using flow cytometry analysis (FACS™ LSRII, BD Biosciences™ Life Sciences™). MoflT1's proliferation analysis of CFSE-stained cells shows distribution of cell populations in 10 generations of cell division. As cells divide, each daughter cell incorporates roughly half of the CFSE dye from its parent allowing cell division to be tracked over time. The higher the generation or of the lower the CFSE fluorescence intensity (i.e., x-axis to the left for the corresponding fluorescence histogram peak). The amplitude or height of each histogram represents the proportion (%) of cells in a given generation.

5 Multi-potency is maintained for hMSCs cultured in BD Falcon Multi-Flasks




hMSCs grown for three consecutive passages in BD Falcon Multi-Flasks retained their multi-potent phenotype and were able to differentiate to adipocytes and osteocytes similar to cells grown in control T-175 flasks. Post-differentiation, cultures were washed with 1X PBS and fixed with 4% paraformaldehyde. Cells cultured for adipogenic differentiation were stained with Oil-Red-O (left panels). Intracellular lipid vacuoles (red) were observed in induced cells and were absent in un-induced control cells. Cells induced for osteogenic differentiation were stained with a substrate for alkaline phosphatase (red) and with a nuclear DAPI stain (blue) (right panel). Alkaline phosphatase positive cells were detected in induced cells but not in un-induced control cells.

3 Immunophenotype of hMSCs Cultured in Media with Serum



Immunophenotypic properties of hMSCs grown in BD Falcon Multi-Flasks for three consecutive passages in serum containing media were nearly identical to cells grown in parallel T-175 control flasks. FACS analysis using BD FACSCanto™ confirmed that hMSCs were negative (filled blue histograms) for surface markers CD14, CD34, CD74, and PL-2 (BD Biosciences™) and positive (filled red histograms) for expression of surface markers CD73, CD90 and CD105 (BD Biosciences). Respective isotype controls for each marker antibody are shown in grey (empty histograms). Cells from three independent flasks of each type were analyzed.

4 Immunophenotype of hMSCs Cultured in Serum-Free Media



Immunophenotypic properties of hMSCs grown in BD Falcon Multi-Flasks for three consecutive passages in serum-free hMSC hMSC SF CCE media were also nearly identical to cells grown in parallel T-175 control flasks. hMSCs were negative (filled blue histograms) for surface markers CD14, CD34, CD74 and PL-2 (BD Biosciences™) and positive (filled red histograms) for expression of surface markers CD73, CD90 and CD105 (BD Biosciences™). Respective isotype controls for each marker antibody are shown in grey (empty histograms). Cells from three independent flasks of each type were analyzed for FACS analysis.

6 Diverse Cell Types Can be Scaled-up Using BD Falcon Multi-Flasks

Cell Type	BD T-175	BD Falcon Multi-Flask (3-Layer)	% Diff.
hMSC-1	2.47 × 10 ⁶	2.46 × 10 ⁶	100%
hMSC-2	3.98 × 10 ⁶	3.93 × 10 ⁶	100%
LuPaq	5.98 × 10 ⁶	5.93 × 10 ⁶	100%
hMSP3	8.98 × 10 ⁶	8.93 × 10 ⁶	100%
hMSP3+325 (PGE2 Inhibitor)	4.98 × 10 ⁶	4.31 × 10 ⁶	86%
CHO-K1	3.98 × 10 ⁶	3.31 × 10 ⁶	83%
hMSC cells from	1.32 × 10 ⁶	1.37 × 10 ⁶	104%
hMSC cells from	3.98 × 10 ⁶	3.93 × 10 ⁶	100%

Equivalent cell yield per unit area was observed with diverse cell types. Data shown is for 3-layer BD Falcon Multi-Flasks. Equivalent cell yield per unit area also was observed for 5-layer BD Falcon Multi-Flasks (data not shown).

For Research Use Only. Not for use in diagnostic or therapeutic procedures. See reader instructions for details. © 2012 BD Biosciences. All rights reserved. BD, the BD logo, and the BD Falcon logo are trademarks of BD Biosciences. All other trademarks are the property of their respective owners. BD Biosciences and Company 2012 001 001