



BD™ CS&T Beads

Catalog No.	Tests
656504	50
656505	150

10/2018

23-14666-03



© 2018 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



Becton, Dickinson and Company
BD Biosciences
2350 Qume Drive
San Jose, CA 95131 USA



Benex Limited
Pottery Road, Dun Laoghaire,
Co. Dublin, Ireland
Tel +353.1.202.5222
Fax +353.1.202.5388

BD Biosciences
European Customer Support
Tel +32.2.400.98.95
Fax +32.2.401.70.94
help.biosciences@europe.bd.com

Australian and New Zealand Distributors:

Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde NSW 2113, Australia

Becton Dickinson Limited
14b George Bourke Drive
Mt Wellington
Auckland 1060, New Zealand

bdbiosciences.com
ClinicalApplications@bd.com

1. INTENDED USE

BD™ CS&T Beads are used on a BD flow cytometer to provide a standardized method to perform quality control of the instrument's optics, electronics, and fluidics, and for adjusting fluorescence compensation. On some BD instruments, BD CS&T Beads are also used for adjusting detector voltages.

2. SUMMARY AND EXPLANATION

BD CS&T Beads are a suspension of fluorospheres with uniform and stable size and fluorescence intensity. The beads are used for instrument quality control (QC) to characterize, track, and report performance measurements of supported flow cytometers. The cytometer's software displays current bead data in plots. Forward scatter (FSC) and side scatter (SSC) identify bead populations based on relative size.

The beads enable the software to measure detector performance and are used to measure the sensitivity of each fluorescence detector. Sensitivity is a measure of the cytometer's ability to resolve dimly stained cells. In addition, the beads are used to optimize the compensation settings each time instrument QC is run.

The software calculates the bright bead median, bright beads %rCV (robust coefficient of variation), and instrument sensitivity for FSC, SSC, and each fluorescence parameter, and compares them to expected values for the bead lot. The rCV measures cytometer alignment.

Daily measurements are automatically entered into Levey-Jennings plots. This allows you to monitor instrument

performance measurements over time and detect potential problems.

For the BD FACSLyric™ flow cytometer, BD CS&T Beads are also used for adjusting detector voltages.

3. PRINCIPLES OF THE PROCEDURE

BD CS&T Beads consist of equal quantities of 3- μm bright, 3- μm mid, and 2- μm dim polystyrene beads. The beads are dyed with fluorochromes. Fluorescence intensity is measured by the cytometer's detectors, processed by the electronics, and displayed and analyzed by the software.

The cytometer's software reads the bead lot file and calculates the target marker position or target fluorescence intensity for the beads. It then calculates the median and %rCV for bright beads.

4. REAGENTS

Reagents provided

BD CS&T Beads are supplied in phosphate buffered saline (PBS) with bovine serum albumin (BSA) and 0.1% sodium azide as follows.

- Two 3-mL vials (50 tests)
- Six 3-mL vials (150 tests)

Each 3-mL vial contains sufficient beads to run 25 tests.

Reagents or materials required but not provided

- Disposable 12 \times 75-mm capped polystyrene test tubes
- Filtered deionized (DI) water, to dilute the beads (BD FACSVia™ system only)
- BD FACSThrow™ sheath fluid (Catalog No. 342003) or equivalent, to dilute the

beads (BD FACSLyric flow cytometer only)

Precautions

- Avoid exposing BD CS&T Beads to direct light.
- Do not run BD CS&T Beads without first diluting them with the proper diluent, as directed in the Procedure section.
- Do not use BD CS&T Beads beyond their expiration date or beyond the day-of-use stability period after dilution. Beads used beyond their stability period begin to lose fluorescence, which can result in failed instrument QC.
- BD CS&T Beads contain sodium azide as a preservative.
- Visit regdocs.bd.com to download the Safety Data Sheet.

Storage and handling

- Store vials at 2°C–8°C and protect from light. Do not use after the expiration date shown on the label.

5. INSTRUMENTS

BD CS&T Beads are for use on the following:

- BD FACSVia flow cytometer
- BD FACSLyric flow cytometer

6. PROCEDURE

Adding or importing bead lot information

Add bead lot ID information by scanning the bead lot file card in this kit.

If you do not have a barcode scanner, import bead lot ID information from the BD Biosciences website.

1. Visit bdbiosciences.com and select **Support** from the menu bar.
The **Services** web page opens.
2. From **Top Support Links** in the right panel, select **Bead Lot Files**: for the appropriate software.
3. Follow the installation instructions on the website to download and import the appropriate bead lot file into the software.
3. Prepare diluted beads according to Table 1 for the system and application you are running.

NOTE Avoid dripping the beads down the side of the tube when diluting them. This can lead to low bead counts during acquisition.

NOTE Do not dilute BD CS&T Beads more than recommended.

4. Vortex the tubes gently before use.

After dilution, the beads are stable for:

- 8 hours at 2°C–25°C on the BD FACSVia system
- 8 hours at 15°C–25°C, or 24 hours at 2°C–8°C on the BD FACSLyric flow cytometer.

NOTE Keep the diluted bead suspension protected from light. Some of the dyes used to manufacture the beads are light sensitive. Fluorescence levels can change if beads are exposed to direct light for longer than 20 minutes.

Preparing a BD CS&T bead suspension

Carefully read the Precautions and Storage and handling statements in the Reagents section.

To prepare the BD CS&T Beads for acquisition:

1. Label a 12 × 75-mm capped polystyrene tube.
2. Thoroughly mix the BD CS&T Beads vial. Invert the vial 10 times, or vortex the vial at medium speed for 5–10 seconds.

Table 1 BD CS&T Beads preparation

For...		Add...				
Task	Cytometer	Diluent	Diluent volume (µL)	Beads (No. of drops)	To tube labeled	How often
Instrument QC (IQC)	BD FACSVia	Filtered DI water	500	2	IQC	<ul style="list-style-type: none"> • Daily • When recommended by BD • To transfer bead lots

Table 1 BD CS&T Beads preparation

For...		Add...				
Task	Cytometer	Diluent	Diluent volume (µL)	Beads (No. of drops)	To tube labeled	How often
Performance QC ^a (PQC)	BD FACSLyric	BD FACFlow sheath fluid	500	2	PQC	Daily
Update reference settings			500	2	Ref	Every 60 days
Bead lot transfer			500	2	Old lot	Before using a new lot
			500	2	New lot	
Characterization QC (CQC)			1,000	4	CQC	<ul style="list-style-type: none"> • Every 6 months • After service or maintenance • When recommended by BD
Lasers setup						

a. Assay and tube settings are automatically updated when running performance QC.

Performing QC on the instrument using BD CS&T Beads

Run instrument QC according to Table 1. See the instrument’s Instructions For Use (IFU) for instructions on installing a bead lot file and performing any of the tasks outlined in Table 1.

7. RESULTS

Reviewing the Instrument QC Report

The Instrument QC Report contains the cytometer serial number, software version, BD CS&T bead lot information, bright bead median, bright beads %rCV, instrument sensitivity, and a pass or fail result for each parameter. A passing result for every parameter is required for instrument QC to pass. A failure for any parameter results in failure of instrument QC. For troubleshooting any QC messages, see the cytometer’s IFU.

When using the BD FACSVia system, we recommend that you visually review the marker positions for scatter and bright bead peaks and adjust the markers to

surround the bead population, as necessary.

8. LIMITATIONS

- BD CS&T Beads are intended for use with supported flow cytometers and their applicable software.
- BD CS&T Beads should not be used to support quantitative fluorescence measurements in a flow cytometer.
- BD CS&T Beads are for instrument QC and setup only.

9. PERFORMANCE CHARACTERISTICS

Performance of the BD CS&T Beads was established for the BD FACSLyric and BD FACSVia flow cytometers by testing at BD Biosciences laboratories in San Jose, CA, USA.

BD FACSLyric Flow Cytometer

Accuracy (BD FACSLyric flow cytometer)

Assay settings and tube settings were determined three times on each of three

BD FACSLyric flow cytometers using one lot of BD CS&T Beads per instrument (total of three lots). For each fluorescence channel, the bright bead median fluorescence intensity (MFI) value (Actual), generated from the assay setup reports, was compared with the lot-specific bright bead MFI value (Target).

Accuracy was calculated as the percent difference between the bright bead MFI values of the Actual and the bright bead MFI values of the Target. See Table 2.

Table 2 Accuracy of cytometer setup using BD CS&T Beads^a

Parameter	Bright bead MFI		% Difference
	Target	Actual	
FSC	17,997	17,998	0.01
SSC	126,535	126,107	-0.34
FITC	5,477	5,479	0.04
PE	12,877	12,873	-0.03
PerCP-Cy5.5	15,574	15,622	0.31
PE-Cy7	11,883	11,867	-0.13
APC	41,354	41,370	0.04
APC-R700 ^b	28,455	28,486	0.11
APC-Cy7	60,190	60,386	0.33
V450 ^a	6,140	6,158	0.29
V500-C ^a	24,442	24,356	-0.35
BV605 ^a	6,457	6,465	0.12
BV711 ^a	38,812	39,014	0.52
BV786 ^a	75,550	75,556	0.01

- a. The data presented are from one run on one instrument with lyse/ no-wash (LNW) tube settings. Results from lyse/wash (LW) tube settings and the other runs and instruments were similar.
- b. BD Horizon™ APC-R700, BD Horizon™ V450, BD Horizon™ V500-C, BD Horizon Brilliant™ Violet 605, BD Horizon Brilliant™ Violet 711, BD Horizon Brilliant™ Violet 786

Reproducibility (BD FACSLyric flow cytometer)

Instrument CQC was run on two BD FACSLyric flow cytometers. Two operators performed two runs of instrument PQC in duplicate on each instrument every day for a period of eight days using one lot of BD CS&T Beads per instrument. Percent coefficient of variation (%CV) of the median MFI values for each channel in high sensitivity and normal modes was used to verify reproducibility. See Table 3.

Table 3 Reproducibility of BD CS&T Beads (Operator/instrument-to-operator/instrument, day-to-day, tube-to-tube)^a

Parameter	%CV (High sensitivity)	%CV (Normal)
FSC	0.49	0.46
SSC	0.28	0.47
FITC	0.24	0.33
PE	0.28	0.27
PerCP-Cy5.5	0.35	0.39
PE-Cy7	0.42	0.52
APC	0.30	0.21
APC-R700	0.37	0.36
APC-Cy7	0.47	0.42
V450	0.26	0.26
V500-C	0.23	0.28
BV605	0.28	0.24
BV711	0.43	0.48
BV786	0.41	0.57

- a. The data presented are from one instrument. Results from the other instrument were similar.

Repeatability (BD FACSLyric flow cytometer)

Instrument CQC and PQC were each performed ten times on two BD FACSLyric flow cytometers using two

lots of BD CS&T Beads to assess run-to-run repeatability. The %CV of the bright beads %rCV (resolution), Br (background), minimum linearity, maximum linearity, Qr (detector

efficiency), and SDen (standard deviation of the electronic noise) were used to verify run-to-run repeatability. See Table 4.

Table 4 Run-to-run repeatability of instrument CQC using BD CS&T Beads^a

Parameter	Bright bead rCV		SD for Br <100	%CV for Br ≥100	SD for linearity minimum (<500)	%CV for linearity maximum	%CV for Qr	%CV for SDen
	%CV of rCV	SD of rCV <2%						
FSC	NA	0.05	NA	NA	NA	NA	NA	NA
SSC	3.27	NA	NA	NA	NA	NA	NA	NA
FITC	NA	0.03	4.45	NA	18.96	0.13	2.74	2.14
PE	NA	0.00	NA	5.57	7.59	0.17	4.72	3.33
PerCP-Cy5.5	1.06	NA	1.34	NA	29.27	0.22	1.84	1.74
PE-Cy7	1.27	NA	0.00	NA	80.02	0.21	2.61	2.47
APC	NA	0.05	1.23	NA	12.75	0.13	3.46	1.91
APC-R700	2.78	NA	9.65	NA	23.88	0.82	4.44	2.69
APC-Cy7	2.97	NA	NA	9.31	28.71	0.18	5.59	1.98
V450	1.76	NA	NA	7.87	45.20	0.10	4.03	2.32
V500-C	1.52	NA	NA	32.91	12.66	0.19	5.75	0.96
BV605	1.62	NA	1.27	NA	41.79	0.19	3.89	1.25
BV711	1.14	NA	0.74	NA	40.23	0.26	4.26	3.02
BV786	1.33	NA	0.00	NA	24.21	0.18	3.96	2.70

a. The data presented are for one lot of BD CS&T Beads run on one instrument. Results for the other lot of beads and instrument for CQC and PQC were similar.

BD FACSVia Flow Cytometer

Reproducibility (BD FACSVia flow cytometer)

Instrument QC was run on two BD FACSVia flow cytometers. Two operators performed two runs of instrument QC in duplicate on each instrument for a period of five days using two lots of BD CS&T Beads.

The coefficient of variation (%CV) of the bright bead median fluorescence intensity (MFI) values and the %CV or standard deviation (SD) of the bright bead robust

coefficient of variation (rCV) values were used to evaluate reproducibility (including instrument-to-instrument, operator-to-operator, day-to-day, run-to-run, and replicate-to-replicate). See Table 5.

Table 5 Reproducibility of BD CS&T Beads (instrument-to-instrument, operator-to-operator, day-to-day, run-to-run, replicate-to-replicate)

Parameter	BD CS&T bead lot	%CV of bright bead MFI	% Bright bead rCV		
			Mean of rCV	%CV of rCV ^a	SD of rCV ^b
FSC-A	Lot A	1.2	2.2	N/A	0.4
	Lot B	1.1	2.2	N/A	0.4
SSC-A	Lot A	1.8	12.3	10.9	N/A
	Lot B	1.8	12.3	11	N/A
FL1-A	Lot A	1.8	2.5	N/A	0.2
	Lot B	1.7	2.5	N/A	0.2
FL2-A	Lot A	4.4	2.2	N/A	0.3
	Lot B	4.3	2.2	N/A	0.2
FL3-A	Lot A	1.2	3.6	N/A	0.2
	Lot B	1.3	3.3	N/A	0.2
FL4-A	Lot A	2.5	4.9	N/A	0.3
	Lot B	2.4	4.6	N/A	0.3

a. %CV was used to evaluate reproducibility when the mean rCV values were greater than 5.

b. SD was used when the mean rCV values were less than 5.

Repeatability (BD FACSVia flow cytometer)
Instrument QC was run on two BD FACSVia flow cytometers. Two operators performed two runs of instrument QC in duplicate on each instrument for a period of five days using two lots of BD CS&T Beads.

The %CV of the bright bead MFI values and the %CV or SD of the %bright bead rCV values were used to evaluate repeatability (replicate-to-replicate). See Table 6.

Table 6 Repeatability of BD CS&T Beads (replicate-to-replicate)

Parameter	BD CS&T bead lot	%CV of bright bead MFI	% Bright bead rCV		
			Mean of rCV	%CV of rCV ^a	SD of rCV ^b
FSC-A	Lot A	0.3	2.2	N/A	0.1
	Lot B	0.3	2.2	N/A	0.1
SSC-A	Lot A	0.4	12.3	1.4	N/A
	Lot B	0.3	12.3	1.8	N/A
FL1-A	Lot A	0.3	2.5	N/A	0.1
	Lot B	0.4	2.5	N/A	0.1

Table 6 Repeatability of BD CS&T Beads (replicate-to-replicate)

Parameter	BD CS&T bead lot	%CV of bright bead MFI	% Bright bead rCV		
			Mean of rCV	%CV of rCV ^a	SD of rCV ^b
FL2-A	Lot A	0.5	2.2	N/A	0.1
	Lot B	0.4	2.2	N/A	0.1
FL3-A	Lot A	0.5	3.6	N/A	0.1
	Lot B	0.4	3.3	N/A	0.1
FL4-A	Lot A	0.4	4.9	N/A	0.1
	Lot B	0.3	4.6	N/A	0.1

a. %CV was used to evaluate repeatability when the mean rCV values were greater than 5.

b. SD was used when the mean rCV values were less than 5.

TROUBLESHOOTING

Problem	Possible Causes	Solution
No beads detected	Beads not mixed prior to diluting, beads are too dilute, there is debris in the bead suspension, incorrect beads were used, beads diluted in wrong diluent, beads exposed to light	<ol style="list-style-type: none"> Vortex the bead vial. Prepare a fresh suspension of beads. Re-run instrument QC.
	Air bubbles in the flow cell or sheath filter	<ol style="list-style-type: none"> For: <ul style="list-style-type: none"> BD FACSVia, perform a backflush or SIP clean. BD FACSLyric, perform a SIT flush. Vortex the tube. Re-run the tube.
	Sheath filter is not filled with fluid	<ul style="list-style-type: none"> For BD FACSVia, perform the two-month maintenance procedure. For BD FACSLyric, purge the sheath filter.

Problem	Possible Causes	Solution
No beads detected	Clogs within the sample path and fluidic lines	<ol style="list-style-type: none"> For: <ul style="list-style-type: none"> BD FACSVia, perform a backflush or SIP clean. BD FACSLyric, perform a SIT flush. Vortex the tube. Re-run the tube.
	Optics are out of alignment	<ul style="list-style-type: none"> Contact BD Biosciences.

Problem	Possible Causes	Solution
Performance check completed with QC messages	Bead gates and markers are not properly adjusted to encapsulate results	For BD FACSVia, review instrument QC results and adjust the CS&T Bead gates for scatter and fluorescence.
	Values for any measurements used to check cytometer performance are not within parameters required for instrument QC to pass (see Reagents section)	Prepare a fresh suspension of beads and re-run instrument QC.
		<ol style="list-style-type: none"> For: <ul style="list-style-type: none"> BD FACSVia, perform the two-month maintenance procedure. BD FACSLytic, perform the monthly cleaning procedure. Re-run the tube.
		Review the instrument QC report to determine whether the specific warnings impact the assay, then continue.
	Contact BD Biosciences.	
Performance check failure	Value(s) for any of the measurements used to check the cytometer performance are not within parameters required for instrument QC to pass (see Section 7)	<ol style="list-style-type: none"> Prepare a fresh suspension of beads. Re-run the performance check.
		Perform the monthly cleaning procedure.
	Improper ratio of 2 μm /3 μm beads due to inadequate mixing of beads	<ol style="list-style-type: none"> Prepare a fresh suspension of beads. Re-run the performance check. If QC fails again, prepare beads from a new vial and be sure to thoroughly vortex the vial prior to use. See Preparing a BD CS&T bead suspension.

For additional troubleshooting assistance, see the cytometer's IFU or contact your local BD Biosciences representative.

WARRANTY

Unless otherwise indicated in any applicable BD general conditions of sale for non-US customers, the following warranty applies to the purchase of these products.

THE PRODUCTS SOLD HEREUNDER ARE WARRANTED ONLY TO CONFORM TO THE QUANTITY AND CONTENTS STATED ON THE LABEL OR IN THE PRODUCT LABELING AT THE TIME OF DELIVERY TO THE CUSTOMER. BD DISCLAIMS HEREBY ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE AND NON-INFRINGEMENT. BD'S SOLE LIABILITY IS LIMITED TO EITHER REPLACEMENT OF THE PRODUCTS OR REFUND OF THE PURCHASE PRICE. BD IS NOT LIABLE FOR PROPERTY DAMAGE OR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING PERSONAL INJURY, OR ECONOMIC LOSS, CAUSED BY THE PRODUCT.