


Sending Results to an LIS Using BD FACSLink™ Software with BD FACSuite™ and BD FACSDuet™ Software

Quick Reference Guide

This quick reference guide (QRG) contains instructions for sending results to a laboratory information system (LIS) using BD FACSLink™ software embedded within BD FACSuite™ software with the BD FACSDuet Sample Preparation System. BD FACSLink software provides an interface to an LIS. It is pre-installed and configured for you either on your server or on a small component that is connected by an Ethernet cable to both the BD FACSDuet™ computer and the computer on which BD FACSuite™ software is installed.


Workflow Overview




Create and add a user-defined (UD) assay

To use BD FACSLink in BD FACSuite software, you must create an assay with the following characteristics:


- Contains a Statistics View with the Include in Auto-Export checkbox selected in the Statistics View Properties dialog.



- The Automatically export results checkbox in the Auto Export Options dialog is selected.



- Multiple Tube Statistics is selected on your worksheet.



- Contains LIS Assay names and Result names.
- Is a shared assay.

Add the UD assay

- 1 Select the **Library** icon.



- 2 Expand **Assays**, select **User-Defined**, then select your assay.

User-Defined Assays	
Name	Modified Date
L&L UD	11/1/2016

- 3 Select the **Send to LIS** tab and click **Edit**.

- 4 Enter a **BD FACSLink Assay Name** and click **Send**.
– Only letters, numbers and underscores are allowed.

- 5 Check the corresponding Statistics box(es) that you would like to import to your LIS.

- 6 Enter the required names in the **BD FACSLink Result Name** column.

Statistics	BD FACSLink Result Name
<input checked="" type="checkbox"/> All Events, Events (Tube_001)	
<input checked="" type="checkbox"/> All Events, % Parent (Tube_001)	
<input checked="" type="checkbox"/> All Events, % Grandparent (Tube_001)	
<input checked="" type="checkbox"/> All Events, % Total (Tube_001)	
<input checked="" type="checkbox"/> All Events, Statistics, FSC-A, Mean (Tube_001)	
<input checked="" type="checkbox"/> All Events, Statistics, SSC-A, Mean (Tube_001)	

Statistics	BD FACSLink Result Name
<input checked="" type="checkbox"/> All Events, Events (Tube_001)	AE
<input checked="" type="checkbox"/> All Events, % Parent (Tube_001)	AEPP
<input checked="" type="checkbox"/> All Events, % Grandparent (Tube_001)	AEPGP
<input checked="" type="checkbox"/> All Events, % Total (Tube_001)	AEPT
<input checked="" type="checkbox"/> All Events, Statistics, FSC-A, Mean (Tube_001)	AEFSCM
<input checked="" type="checkbox"/> All Events, Statistics, SSC-A, Mean (Tube_001)	AESSCM

A confirmation dialog will be displayed.

Publish the UD FACSuite assay to FACSDuet software

When you have a UD assay in BD FACSuite, publish the assay to BD FACSDuet software to allow for a preparation method to be added.

In BD FACSuite software you must create an assay with the following characteristics:

- Contains an LIS assay name and result name.
- Is a shared assay.

- 1 Click the **Library** icon to go to your library.



- 2 Expand **Assays**, select **User-Defined**, then select your assay.

Library		User-Defined Assays	
		Name	Modified Date
Assays			
BD			
User-Defined		L&L UD	11/1/2016
Beads and Reagents			
Keywords			
Labels			
Tube Settings			

- 3 Right-click on an assay name and select **FACSDuet Publish**.

User-Defined Assays				
Name	Modified Date	Author		
2222 UD	2/4/2019	Admin User		Y
City Of Hope Demo UD	1/31/2019	Admin User		Y
Lymph UD	2/11/2019	Admin User		Y
NK Activated T Cells UD	2/11/2019	Admin User		Y
T Cell Activation Panel UD	2/11/2019	Admin User		Y
Z Test UD		Admin User		Y

Rename

FACSDuet Publish

FACSDuet Export

Note: A message will appear only if the assay cannot be published.


Import the UD FACSuite assay to FACSDuet software

After you have published an UD FACSuite assay to BD FACSDuet software, you must also import the assay.


- 1 In BD FACSDuet software, select the **Library** icon.




- 2 Under the Assays tab, tap **Import**.




- 3 A dialog will appear. Click **FACSLink/Import** from middleware.




- 4 A dialog will appear. Select the assay you want to import.



- 5 Select the assay you imported, then tap the **Edit** icon and continue creating the assay.





- 6 Ensure that each tube has the appropriate preparation method assigned. Select a tube and tap **Options > Choose Method > Apply**. When the assay is complete, tap **Save**.



Send results to the LIS

- 1 After acquisition, analyze the entry and mark the entry **Approved** in either the **Entry Details** panel or the **Worklist Entries** panel. The results are automatically sent to the LIS upon approval of an entry.



- 2 Select **Tools > BD FACSLink Test Results History** to view information about when the results for all worklists were transferred.

The screenshot shows the 'BD FACSLink Test Results History' window. It contains a table with the following data:


Status	Worklist Name	Sample ID	Task Name	Approved Time	Delivered Time
Sent	Worklist_001	B123	3/16+56/45/19 + Tr	10/21/2016 11:51:2:	10/21/2016 11:51:23 AM
Sent	Worklist_001	C123	3/16+56/45/19 + Tr	10/21/2016 11:51:2:	10/21/2016 11:51:22 AM
Sent	Worklist_001	B123	4 Color TBNK + Truc	10/21/2016 11:51:1!	10/21/2016 11:51:20 AM
Sent	Worklist_001	B123	4 Color TBNK + Truc	10/21/2016 11:46:3!	10/21/2016 11:46:31 AM
Sent	Worklist_002	123	6 Color TBNK + Truc	09/13/2016 12:29:5:	09/13/2016 12:29:54 PM

At the bottom of the window, there are three buttons: 'Refresh', 'Delete', and 'Close'.

Test the connections

Test the connection to BD FACSLink software

- 1 Start BD FACSuite Clinical software.
- 2 Select **Tools > Preferences**, then select the **Worklist** tab.
- 3 Verify the following information, then click **Test Connection**.



BD FACSLink Connection

Cytometer Serial Number: Z6591800011

Username: IIM_ADMIN


Password:

TCP/IP Address: http :// 172.16.10.1 Port: 10005 **Test Connection**

Client Certificate: Browse...

Note: Only users with Administrator privileges can change connection settings.


- 4 Click **OK**.
- Check the Troubleshooting table at the end of this QRG for resolving connection problems.



System Status			
Cytometer	Universal Loader	Workstation	BD FACSLink
BD FACSLyric (SN:12345)	Installed	1381KC2	<input checked="" type="checkbox"/> Connected ←


Test the connection to BD FACSDuet software

- 1 Start BD FACSDuet software.
- 2 Select **Settings > System**.
- 3 Verify the following information, then tap **Test Connection**.



Note: Only users with Administrator privileges can change connection settings.



Note: Consult the Troubleshooting table for help with resolving connection problems.



Troubleshooting Connection Issues

As a workaround until you resolve a connection issue with BD FACSLink software, you can export your worklists and import them into BD FACSuite software. See Exporting Worklists in the BD FACSDuet user guide.

Note: Only users with Administrator privileges can change connection settings.

Observation	Possible causes	Recommended solutions
Unable to send orders down or results up to the LIS	Connection error	<ol style="list-style-type: none"> 1. Log into Data Innovations Instrument Manager™ as an Administrator. 2. Select System > Status. 3. Check if the Status column shows Off for either FACSDuet, FACSLyric, or both. <p>If the Status is:</p> <p>Off → </p>  <p>If the Status shows Connecting and then turns off, read the error message and select either Start Selected Connections, Stop Selected Connections, or Resend Messages.</p>
BD FACSLink software is not connected	<p>If BD FACSLink software connects through an Ethernet cable, the cable might be disconnected.</p> <p>The BD FACSLink connection settings are invalid.</p>	<ul style="list-style-type: none"> • Check the Ethernet cable connection on the back of the BD FACSDuet computer and on the BD FACSLink computer. • Check the Ethernet cable connection on the BD FACSLyric computer and on the BD FACSLink computer. <p>After the BD FACSDuet system is installed, it should not be necessary to change the BD FACSLink connection settings. If you suspect that the connection settings need to be changed, contact BD Technical Support. In preparation, do the following steps:</p> <ol style="list-style-type: none"> 1. With Administrator privileges, tap Settings > System and check the connection settings. 2. Ensure that there are no obvious issues that you can easily correct, such as empty fields or wrong instrument IDs. 3. Make a note of the settings so that you can provide this information to BD Technical Support.

This material is for training purposes.

BD FACSuite Clinical Software and BD FACSDuet Software are CE marked For In Vitro Diagnostic Use.

BD FACSLink software is CE marked For In Vitro Diagnostic Use.

BD FACSuite software is for Research Use Only. Not for use in diagnostic or therapeutic procedures.

23-21384-00

Becton, Dickinson and Company, BD Biosciences, 2350 Qume Drive, San Jose, CA 95131 USA

bdbiosciences.com

BD, the BD Logo, FACSDuet, FACSLink, FACSLyric, and FACSuite are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2019 BD. All rights reserved.

