

Flow Cytometry Panel Design Journey

We're with you every step of the way

Waterfall of Discovery

There are hundreds of research applications for flow cytometry.

Define your experimental hypothesis

Step 1: Defining your experimental hypothesis is the first step in panel design. Start with identifying:

- The biological information you are trying to achieve
- The population(s) of cells you wish to interrogate
- Whether targets are found on the cell surface or intracellularly

Did you know?
BD launched the world's 1st commercially available flow cytometer in 1974 - kicking off a legacy of research innovation.

Viability Volcano

Beware of hitchhiking dead cells because they can skew your data. Be sure to identify them properly with a viability dye.

Know your flow cytometer

Step 3: Knowing your instrument is essential. Understanding your instrument's configuration will let you know how many markers and which fluorochromes your instrument can detect.

Elements to consider include:

- Laser wavelength for excitation
- Number of detectors off each laser
- Filters available to detect the fluorochromes

Did you know?
BD offers over 40 fluorochromes and more than 1,500 specificities to provide researchers flexibility in panel design.

Marker selection

Step 2: During the second step of the panel design process, you will need to identify which and how many markers you need to identify the population of interest.

Pay attention to:

- Marker expression levels
 - Primary antigen: Expressed at high density, often defining lineages
 - Secondary antigen: Often expressed over a continuum
 - Tertiary antigen: Critical markers expressed at low density
- Marker coexpression, especially of dim markers
- The gating strategy needed to identify the population(s) of cells you wish to interrogate

Fluorochrome assignment

Step 4: Carefully select fluorochromes to resolve markers at all expression levels and minimize spectral overlap. Consider using tools like a fluorochrome resolution ranking and a spectrum viewer to help assess:

- Cross laser excitation
- Fluorochrome spillover

Remember to pair bright fluorochromes with low expressing antigens and dim fluorochromes with high expressors. Keep in mind that spread only impacts the resolution of coexpressed markers.

Lake of Permeabilization

For smooth sailing, select the appropriate buffer system for your antibodies.

Did you know?
Since the introduction of revolutionary Sirigen polymer dye technology, BD has launched 17 innovative dyes including BD Horizon Brilliant™ UV and BD Horizon Brilliant Violet™ Dyes.

TIP: For help choosing fluorochromes, check out our BD Spectrum Viewer: bdbiosciences.com/spectrum-viewer

Spectral Mountain Passage

Map out the peaks and valleys of the dyes' spectra to minimize spillover and spread.

Review panel

Step 5: Review your panel design and begin ordering your reagents.

Remember to titrate your mass size reagents and optimize your staining protocol. Include proper controls for compensation, FMO and biological controls to help ensure optimal panel performance.

Ready to begin workflow

Check out our online flow cytometry panel design educational resources

-  Videos
-  Webinars
-  Panel design tools



bdbiosciences.com/panel-design

To learn more about flow cytometry panel design resources or for support, please contact your BD Sales Representative.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.
23-23036-00

BD Life Sciences, San Jose, CA, 95131, USA

bdbiosciences.com



BD, the BD Logo, Horizon Brilliant and Horizon Brilliant Violet are trademarks of Becton, Dickinson and Company or its affiliates.
© 2020 BD. All rights reserved.