

BD Phosflow Human CrkL (pY207) Monoclonal Antibodies

Features

Useful for the identification of cells that express the BCR-ABL fusion protein

Suitable marker for the study of CrkL signal transduction

Available as purified and in a wide selection of fluorochrome-conjugated formats including Alexa Fluor® 488, Alexa Fluor® 647, and PE

BD Biosciences now offers purified and fluorochrome-conjugated antibodies for the study of phosphorylated Crk-Like (CrkL), a key adaptor molecule in BCR-ABL signaling.

Anti-human CrkL (pY207) antibody (clone K30-391.50.80)

The new BD™ Phosflow human CrkL (pY207) antibody (clone K30-391.50.80) reacts with CrkL when phosphorylated at tyrosine 207. It does not react with the unphosphorylated protein. CrkL is a 39-kDa adaptor protein that is preferentially expressed in hematopoietic cells. CrkL contains SH2 and SH3 binding domains that allow it to interact with a variety of effector proteins including paxillin, p130Cas, c-Cbl, c-Abl, and C3G.

Available conjugates include Alexa Fluor® 488, Alexa Fluor® 647, and PE formats to enable maximum flexibility for design of multicolor panels in combination with any of our family of BD FACST™ brand flow cytometers. Antibodies were tested in human model systems but are likely to cross-react with mouse and rat model systems due to predicted sequence identity.

Discovery of Crk proteins


Crk adaptor proteins were discovered in the 1980s as a novel retroviral gene product v-Crk. CrkL is a cellular homolog to v-Crk and contains Src Homology 2 (SH2) and Src Homology 3 (SH3) domains separated by flexible linker sequences. SH2 and SH3 domains bind to many proteins and aid in the formation of molecular complexes. SH2 domains typically bind to proteins that contain a phosphorylated tyrosine residue followed by specific residues.^{1,2}

CrkL proteins are important for many biological processes including cell adhesion, cell proliferation, cell migration, apoptosis, and regulation of gene expression.¹

The role of CrkL in BCR-ABL signaling

The human gene encoding the CrkL protein was discovered through its localization on chromosome 22 near the breakpoint cluster region (BCR). CrkL forms complexes with BCR-ABL. CrkL interacts with the proline-rich region of the C-terminal of the Abl kinase.²

Mutations of BCR-ABL in the proline rich C-terminal, as well as the Grb2 and CBL binding domains, prevent the induction of leukemia in mice by BCR-ABL. These mutations of the BCR-ABL pathway severely compromise the ability of BCR-ABL to activate PI3 kinase and MAP kinase pathways. Interestingly, Stat5 activation is maintained.³



Overlay Gate	Cell Population	Median
Gate 1	CD3 ⁺	8.98
Gate 2	CD34 ⁺	93.06

Figure 1. Peripheral blood mononuclear cells from a qPCR BCR-ABL-positive sample were fixed, permeabilized, and stained with PerCP anti-CD3, PE anti-CD34, and Alexa Fluor® 647 anti-phospho-CrkL (Y207) clone K30-391.50.80 (Cat. No. 560790). Phospho-CrkL staining was assessed on CD3-positive and CD34-positive cells. The CD34-positive cells have 10-fold greater phospho-CrkL than the CD3 cells (negative control) from the same sample. Data courtesy of Dr. Brian Druker.

Visit bdbiosciences.com/phosflow for more information.

BD™ Phosflow Human CrkL (pY207) Monoclonal Antibodies

References

1. Birge RB, Kalodimos C, Inagaki F, Tanaka S. Crk and CrkL adaptor proteins: networks for physiological and pathological signaling. *Cell Commun Signal*. 2009;7:13-36.
2. Feller SM. Crk family adaptors-signalling complex formation and biological roles. *Oncogene*. 2001;20:6348-6371.
3. Johnson KJ, Griswold IJ, O'Hare T, et al. A BCR-ABL mutant lacking direct binding sites for the GRB2, CBL and CRKL adaptor proteins fails to induce leukemia in mice. *PLoS One*. 2009;4:e7439.

Ordering Information

Description	React	Clone	Isotype	Format	Size	Cat. No.
CrkL (pY207)	Hu	K30-391.11.30	Mouse IgG _{2a} , κ	Purified	0.1 mg	558386
		K30-391.50.80	Mouse IgG _{2a} , κ	Alexa Fluor® 488	50 tests	560789
				Alexa Fluor® 647	50 tests	560790
				PE	50 tests	560788

Related Products

Description	React	Clone	Isotype	Format	Size	Cat. No.
BCR-ABL Protein Kit				BD™ CBA Kit	50 tests	643939
BCR (pY177)	Ms	J52-309	Mouse IgG _{2b} , κ	Purified	0.1 mg	558248
Stat5 (PY694)	Hu, Ms	47	Mouse IgG ₁ , κ	Alexa Fluor® 488	50 tests	612598
				Alexa Fluor® 647	50 tests	612599
				Pacific Blue™	50 tests	560311
				PE	50 tests	612567
				PE-Cy™7	50 tests	560117
				PerCP-Cy™5,5	50 tests	560118
BD Phosflow Perm Buffer III					125 mL	558050
BD Phosflow Perm Buffer IV (10X)					50 mL	560746

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

Alexa Fluor® is a registered trademark and Pacific Blue is a trademark of Molecular Probes, Inc. The Alexa Fluor® and Pacific Blue™ dye antibody conjugates in this product are sold under license from Molecular Probes, Inc., for research use only, excluding use in combination with microarrays, or as analyte specific reagents. The Alexa Fluor® dyes (except for Alexa Fluor® 430) and Pacific Blue™ dye are covered by pending and issued patents.

Cy™ is a trademark of Amersham Biosciences Corp. Cy™ dyes are subject to proprietary rights of Amersham Biosciences Corp. and Carnegie Mellon University and are made and sold under license from Amersham Biosciences Corp. only for research and in vitro diagnostic use. Any other use requires a commercial sublicense from Amersham Biosciences Corp., 800 Centennial Avenue, Piscataway, NJ 08855-1327, USA.

BD, BD Logo and all other trademarks are property of Becton, Dickinson and Company. © 2010 BD

23-11332-00



BD Biosciences

2350 Qume Drive
San Jose, CA 95131
US Orders: 877.232.8995
answers@bd.com
bdbiosciences.com