

Job Aid

BD FACSDiscover™ Family: Managing data

This job aid contains instructions for how to export, import, and delete data for the BD FACSDiscover™ Family. You can export a Chorus Experiment File (CEF), flow cytometry data (FCS), imaging data (CVW), and statistics (CSV).

NOTE The exported FCS data and statistics will reflect the unmixing matrix (recording specific or latest) that is selected in the experiment at the time of export.

For additional information, see the *BD FACSDiscover™ A8 Cell Analyzer or BD FACSDiscover™ S8 Cell Sorter user's guides*.




Before you begin

Power on the instrument workstation and log into BD FACSDiscover™ Software.

Managing data from the Experiments page

Exporting data

1. Navigate to the Experiments page.
2. Select the checkboxes for the experiments you need to export.
3. Export the data.
 - a. Click **Export**.
 - b. Click **Change** to select a location on the workstation to save the files.



Exporting data, continued

c. Select the files and data formats to be exported.

<input checked="" type="checkbox"/>	All				<input checked="" type="checkbox"/> Select All Chorus Experiment Files (1)		
<input checked="" type="checkbox"/>	▼ Plate_001 (1)						
	Name	Events	Created By	Size	<input checked="" type="checkbox"/> FCS	<input checked="" type="checkbox"/> Statistics	<input checked="" type="checkbox"/> Images
<input checked="" type="checkbox"/>	A01-Sample_001	92,000	admin user	8 GB	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	▼ Plate_002 (1)						
	Name	Events	Created By	Size	<input checked="" type="checkbox"/> FCS	<input checked="" type="checkbox"/> Statistics	<input checked="" type="checkbox"/> Images
<input checked="" type="checkbox"/>	A01-Sample_001	141,000	admin user	12 GB	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	▼ Rack_002 (1)						
	Name	Events	Created By	Size	<input checked="" type="checkbox"/> FCS	<input checked="" type="checkbox"/> Statistics	<input checked="" type="checkbox"/> Images
<input checked="" type="checkbox"/>	A1-Sample_001	17,000	admin user	1 GB	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>

TIP To reduce the size of the exportable FCS file, select a size from the **FCS Parameter Selection** menu.

☒ Select All Chorus Experiment Files (1)

☒ Include Chorus Experiment File (1)

FCS Parameter Selection

Recommended (Approx. 30% reduction)
Recommended (Approx. 30% reduction)
Small (Approx. 60% reduction)
All

d. Click **Export**.

When data export is complete, close the dialog.

6 item(s) selected

Close

Export

Deleting data

1. Navigate to the Experiments page.
2. Select the checkboxes for the experiments where you need to delete data.
3. Delete the data.

a. Click **Delete Data**.

b. Select the files to be deleted.

c. Click **Delete Selected**.

Import

Export (12)

Delete Data (5)

Delete Experiment (1)

Delete Selected

11 item(s) selected

11 item(s) selected

4. Close the dialog, if needed.

Deleting experiments


- 1. Navigate to the Experiments page.
- 2. Select the checkboxes for the experiments you need to delete.

Cytometer		Name	Mode	Tube(s)	Carrier(s)
	<input checked="" type="checkbox"/>	5 Color PBMC MM	Imaging	2	1

- 3. Click **Delete Experiment**.




- 4. Click **Yes** in the dialog.
- NOTE** Deleting an experiment that is used as a template will also delete the template from the New Experiment dialog.




Importing experiments

- 1. Navigate to the Experiments page.
- 2. Click **Import**.
- 3. Browse to the folder on the workstation that contains the Chorus Experiment File (CEF) of interest.
- 4. Select the CEF file and click **Open**.



- 5. (Optional) Clear the **Include Data** checkbox.
- NOTE** If the CEF file is in a folder that also contains FCS data, the checkbox is selected by default.



- 6. Click **Import**.
- If a copy of the experiment exists, the name of the imported experiment is modified automatically.

Users				
	<input type="checkbox"/>	6 color sort_AK 1	6	
	<input type="checkbox"/>	6 color sort_AK	6	


Extracting images

Extracting TIF files from CVW files

1. Open the BD CellView™ Image Extractor.

The image extractor was installed on the workstation desktop during instrument installation. It is also available for download through the BD FACSCorus™ Software Help menu.

TIP Use a standalone computer when extracting images for better computer performance.



2. Click **Select Folder**.

3. Browse to the folder on the workstation that contains the CVW files of interest and select it, then click **Select Folder**.


TIP Selecting the experiment folder is sufficient to extract images from all the data files within the experiment.

4. Click the checkboxes to select the CVW files to be converted.

5. Click **Extract Images**.

6. (Optional) Click **Change Location** to select another folder containing CVW files.

7. Close the BD CellView™ Image Extractor.



This material is for training purposes.
For Research Use Only. Not for use in diagnostic or therapeutic procedures.

BD Life Sciences, Milpitas, California, 95035, USA

bdbiosciences.com

BD, the BD Logo, BD CellView, BD FACSCorus, BD FACSD Discover and BD SpectralFX are trademarks of Becton, Dickinson and Company or its affiliates. ©2025 BD. All rights reserved. NPM-6380 (v2.0) 1125

