



Analyze extracellular vesicles (EVs) on the BD FACSymphony™ A1 Cell Analyzer


Standardized workflow for consistent results



Fine tune your small particle detector (SPD)

STEP 1: Prepare your instrument using standard small particle beads

- Prepare SPD bead mixture
- Align SPD using Picomotor Application (if needed)
- Adjust SP SSC voltage to place bead peaks on target



Ensure reproducible low background noise

STEP 2: Perform start-up procedure to reduce background noise

- Run BD® Detergent Solution Concentrate
- Run sample diluent and check background noise


3.1 Buffer only*	3.5 Single-stained controls*
3.2 Buffer with reagents* 3.6	Procedural controls*
3.3 Unstained controls*	3.7 Serial dilution*
3.4 Isotype controls	3.8 Detergent-treated EV samples

Welsh JA, et al. Journal of Extracellular Vesicles. 2020;9(1):1713526.*

Acquire your EV samples and explore the power of “OR” thresholding

STEP 3: Apply proper controls, calibration and thresholding strategy

- Use *MiFlowCyt-EV* framework for publication to guide you on assay controls and SP SSC and fluorescence calibration
- Use “OR” thresholding for SP SSC and relevant fluorescent parameters for maximum information
- Consider running serial dilution series to check for absence of swarming



Analyze your EV data

STEP 4: Take advantage of the analysis power of BD FACSDiva™ and FlowJo™ Software

- Output files from recommended calibration software compatible with BD FACSDiva™ and FlowJo™ Software for further analysis

Still need help?

- Follow troubleshooting guidelines in the User's Guide
- Contact your BD support team for assistance

To learn more about the BD FACSymphony™ A1 Cell Analyzer or the BD® Small Particle Detector, contact your BD sales representative.

*User should follow routine maintenance process and run CS&T before SPD procedure
Class 1 Laser Product.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

BD Life Sciences, San Jose, CA 95131, USA

bdbiosciences.com

