

BD CD10 (HI10a)

Form	Catalog number	Form	Catalog number
FITC	340924	APC	340922
PE	340920	APC-H7	655426
PerCP-Cy5.5	663194		
PE-Cy7	341102		

BD Horizon™ CD10 (HI10a)

Form	Catalog number
APC-R700	659111
BV421	659449
BV605	663986
BV711	664520

Description

Specificity

The CD10 antibody recognizes a 100-kilodalton (kDa) type II transmembrane, glycosylated, zinc-containing metalloprotease.^{1,2} The CD10 antigen is also known as common acute lymphoblastic leukemia antigen (CALLA), neutral endopeptidase (NEP), gp100, and enkephalinase.³

Antigen distribution

The CD10 antigen is found on lymphocytes from samples with acute B-lymphoid leukemia.⁴ The CD10 antigen is also present on a wide variety of normal and neoplastic cell types including renal epithelium, fibroblasts, granulocytes, germinal center B lymphocytes,⁵ neutrophils,⁶⁻⁸ some T-cell leukemias,⁹ and some lymphoma, melanoma, and glioma cell lines.³

The CD10 antigen cleaves a number of biologically active peptides,¹⁰ including fMLP, and may modulate the chemotactic activity of fMLP towards neutrophils.¹¹ Inhibition of the CD10 antigen promotes B-cell maturation,¹² suggesting that it plays a role in B-cell development.

Clone

The CD10 antibody, clone HI10a,² is derived from the hybridization of P3-63-Ag8.653 mouse myeloma cells with spleen cells isolated from BALB/c mice immunized with blasts from a patient with acute CALLA leukemia.

Composition

The CD10 antibody is composed of mouse IgG₁ heavy chains and kappa light chains.

Product configuration

The following are supplied in buffer containing a stabilizer and a preservative.

Analyte Specific Reagent. Analytical and performance characteristics are not established.

Form	Number of tests	Volume per test (μL)	Amount provided (μg)	Total volume (mL)	Concentration (μg/mL)	Stabilizer	Preservative
FITC	50	20	12.5	1	12.5	Gelatin	0.1% Sodium azide
PE	50	20	6	1	6	Gelatin	0.1% Sodium azide
PerCP-Cy5.5	100	5	12.5	0.5	25	Gelatin	0.1% Sodium azide
PE-Cy7	100	5	12.5	0.5	25	Gelatin	0.1% Sodium azide
APC	100	5	12.5	0.5	25	Gelatin	0.1% Sodium azide
APC-H7	100	5	25	0.5	50	BSA	0.0026% CMIT/MIT (3:1)
APC-R700 ^a	100	5	6.25	0.5	12.5	BSA	0.0028% CMIT/MIT (3:1)
BV421 ^a	100	5	3.15	0.5	6.3	BSA	0.009% MIT
BV605 ^a	100	5	25	0.5	50	BSA	0.09% Sodium azide
BV711 ^a	100	5	25	0.5	50	BSA	0.09% Sodium azide
^a BD Horizon™ APC-R700, BD Horizon Brilliant™ Violet 421, BD Horizon Brilliant™ Violet 605, BD Horizon Brilliant™ Violet 711							

CAUTION Some PE-Cy7, APC-H7, and APC-R700 conjugates show changes in their emission spectra with prolonged exposure to paraformaldehyde or light. For overnight storage of stained cells, wash and resuspend in buffer without paraformaldehyde after 1 hour of fixation.

CAUTION Prolonged exposure of cells to paraformaldehyde can lead to increased autofluorescence in the violet channels. For overnight storage of stained cells, wash and resuspend in buffer without paraformaldehyde after 1 hour of fixation.

CAUTION If you choose to combine BD Horizon Brilliant™ reagents in a multicolor staining cocktail, dyes might bind to one another without the use of a buffering solution, such as BD Horizon™ Brilliant Stain Buffer.

Purity

FITC: ≤5% free fluorophore, based on size-exclusion chromatography (SEC)

PE, PerCP-Cy5.5, PE-Cy7, APC, APC-H7, APC-R700: ≤20% free fluorophore at bottling, as measured by SEC

BV421, BV605, BV711: ≤25% free fluorophore, as measured by ion-exchange chromatography (IEC)


Handling and Storage

Store vials at 2–8 °C. Conjugated forms should not be frozen. Protect from exposure to light. Each reagent is stable until the expiration date shown on the bottle label when stored as directed.

Warning

All biological specimens and materials coming in contact with them are considered biohazards. Handle as if capable of transmitting infection^{13,14} and dispose of with proper precautions in accordance with federal, state, and local regulations. Never pipette by mouth. Wear suitable protective clothing and gloves.

The APC-H7 contains 0.0026% and APC-R700 contains 0.0028% of a mixture of 5-Chloro-2-methyl-4-isothiazolin-3-one and 2-Methyl-2H -isothiazol-3-one [CMIT/MIT (3:1)], CAS number 55965-84-9 and BV421 contains 0.009% 2-Methyl-4-isothiazolin-3-one, CAS number 2682-20-4. These reagents are classified as hazardous according to the Globally Harmonized System of Classification and Labelling of Chemicals (GHS).

	Warning
	H317: May cause an allergic skin reaction. H412: Harmful to aquatic life with long lasting effects.
Prevention	P261: Avoid breathing dust/fume/gas/mist/vapors/spray. P272: Contaminated work clothing should not be allowed out of the workplace. P273: Avoid release to the environment. P280: Wear protective gloves/protective clothing/eye protection/face protection.
Response	P302+P352: IF ON SKIN: Wash with plenty of soap and water. P333+P313: If skin irritation or rash occurs: Get medical advice/attention. P321: Specific treatment (see supplemental first aid instructions on this label). P362+P364: Take off contaminated clothing and wash it before reuse.
Disposal	P501: Dispose of contents/ container to an approved facility in accordance with local, regional, national and international regulations.

Go to regdocs.bd.com to download the Safety Data Sheet.

Characterization

To ensure consistently high-quality reagents, each lot of antibody is tested for conformance with characteristics of a standard reagent.

Warranty

Unless otherwise indicated in any applicable BD general conditions of sale for non-US customers, the following warranty applies to the purchase of these products.

THE PRODUCTS SOLD HEREUNDER ARE WARRANTED ONLY TO CONFORM TO THE QUANTITY AND CONTENTS STATED ON THE LABEL OR IN THE PRODUCT LABELING AT THE TIME OF DELIVERY TO THE CUSTOMER. BD DISCLAIMS HEREBY ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY AND FITNESS FOR ANY PARTICULAR PURPOSE AND NONINFRINGEMENT. BD'S SOLE LIABILITY IS LIMITED TO EITHER REPLACEMENT OF THE PRODUCTS OR REFUND OF THE PURCHASE PRICE. BD IS NOT LIABLE FOR PROPERTY DAMAGE OR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING PERSONAL INJURY, OR ECONOMIC LOSS, CAUSED BY THE PRODUCT.

References

1. Greaves MF. Monoclonal antibodies as probes for leukemic heterogeneity and hematopoietic differentiation. In: Knapp W, ed. *Leukemia Markers*. New York, NY: Academic Press; 1981:19.
2. Zola H. CD10 Workshop Panel report. In: Schlossman SF, Boumsell L, Gilks W, et al, eds. *Leucocyte Typing V: White Cell Differentiation Antigens*. Vol 1. New York, NY: Oxford University Press; 1995:505-507.
3. Letarte M, Vera S, Tran R, et al. Common acute lymphocytic leukemia antigen is identical to neutral endopeptidase. *J Exp Med*. 1988;168:1247-1253.

4. LeBien TW, McCormack RT. The common acute lymphoblastic leukemia antigen (CD10): emancipation from a functional enigma. *Blood*. 1989;73:625-635.
5. Caligaris-Cappio F, Riva M, Tesio L, Schena M, Gaidano G, Bergui L. Human normal CD5⁺ B lymphocytes can be induced to differentiate to CD5⁻ B lymphocytes with germinal center cell features. *Blood*. 1989;73:1259-1263.
6. Connelly JC, Chambless R, Holiday D, Chittenden K, Johnson AR. Up-regulation of neutral endopeptidase (CALLA) in human neutrophils by granulocyte-macrophage colony-stimulating factor. *J Leukoc Biol*. 1993;53:685-690.
7. Connelly JC, Skidgel RA, Schulz WW, Johnson AR, Erdos EG. Neutral endopeptidase 24.11 in human neutrophils: cleavage of chemotactic peptide. *Proc Natl Acad Sci USA*. 1985;82:8737-8741.
8. Erdös EG, Wagner B, Harbury CB, Painter RG, Skidgel RA, Fa XG. Down-regulation and inactivation of neutral endopeptidase 24.11 (enkephalinase) in human neutrophils. *J Biol Chem*. 1989;264:14519-14523.
9. Consolini R, Legitimo A, Rondelli R, et al. Clinical relevance of CD10 expression in childhood ALL. *Haematologica*. 1998;83:967-973.
10. Erdös EG, Skidgel RA. Neutral endopeptidase 24.11 (enkephalinase) and related regulators of peptide hormones. *FASEB J*. 1989;3:145-151.
11. Hofman P, Selva E, Le Negrate G, et al. CD10 inhibitors increase f-Met-Leu-Phe-induced neutrophil transmigration. *J Leukoc Biol*. 1998;63:312-320.
12. Salles G, Rodewald HR, Chin BS, Reinherz EL, Shipp MA. Inhibition of CD10/neutral endopeptidase 24.11 promotes B-cell reconstitution and maturation in vivo. *Proc Natl Acad Sci USA*. 1993;90:7618-7622.
13. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI document M29-A4.
14. Centers for Disease Control and Prevention. 2007 Guideline for Isolation Precautions: Preventing Transmission of Infectious Agents in Healthcare Settings. <https://www.cdc.gov/infectioncontrol/guidelines/isolation/index.html>. Accessed March 12, 2019.

Patents and Trademarks

For U.S. patents that may apply, see [bd.com/patents](https://www.bd.com/patents).

BD, the BD Logo, BD Horizon Brilliant and Horizon are trademarks of Becton, Dickinson and Company or its affiliates. All other trademarks are the property of their respective owners. © 2023 BD. All rights reserved.

Cy™ is a trademark of GE Healthcare. This product is subject to proprietary rights of GE Healthcare and Carnegie Mellon University, and is made and sold under license from GE Healthcare. This product is licensed for sale only for in vitro diagnostics. It is not licensed for any other use. If you require any additional license to use this product and do not have one, return this material, unopened, to BD Biosciences, 155 North McCarthy Boulevard, Milpitas, California 95035, USA, and any money paid for the material will be refunded.

Contact Information

Becton, Dickinson and Company

BD Biosciences

155 North McCarthy Boulevard

Milpitas, California 95035 USA

bdbiosciences.com

ClinicalApplications@bd.com



eifu.bd.com