

**BD OptEIA™**

# Human IL-2 ELISA Kit II

## Instruction Manual

Catalog No. 550611



*BD flow cytometers are Class 1 Laser Products.*

*For Research Use Only. Not for use in diagnostic or therapeutic procedures.*

*© 2018, Becton, Dickinson and Company. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without prior written permission from BD Biosciences.*

*The information in this guide is subject to change without notice. BD Biosciences reserves the right to change its products and services at any time to incorporate the latest technological developments. Although this guide has been prepared with every precaution to ensure accuracy, BD Biosciences assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information.*

*BD Biosciences welcomes customer input on corrections and suggestions for improvement.*

*© 2018 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.*

# Table of Contents

Introduction.....	4
Principle of the Test.....	4
Reagents Provided.....	5
Materials Required but not Provided.....	5
Storage Information.....	6
Warnings and Precautions.....	6
Specimen Collection and Handling.....	8
Reagent Preparation.....	8
Assay Procedure.....	10
Assay Procedure Summary.....	11
Calculation of Results.....	11
Typical Data.....	12
Limitations of the Procedure.....	12
Performance.....	13
Limit of Detection.....	13
Recovery.....	13
Linearity.....	13
Specificity.....	14
Precision.....	15
Intra-assay.....	15
Inter-assay.....	15
Standardization.....	15
Experimental Results.....	15
Serum.....	15
Plasma.....	16
Cell Culture Supernatants.....	16
Troubleshooting.....	17
References.....	18
Plate Templates.....	19

# Introduction

Interleukin-2 (IL-2), originally called T cell growth factor (TCGF), is a multi-functional cytokine that can stimulate the differentiation and proliferation of T lymphocytes and other cell types including B lymphocytes, NK cells, LAK cells, and monocytes/macrophages.<sup>1-3</sup> IL-2 exerts its biological effects by binding to specific receptors expressed by various target cells. Human IL-2 is a 15 kDa protein containing 133 amino acid residues.

The BD OptEIA™ ELISA Kit II format was developed for superior accuracy with serum and plasma specimens. The data that demonstrates this enhancement can be located in the Performance “Recovery” and “Linearity” sections.

The BD OptEIA™ Human IL-2 ELISA Kit II is for the quantitative determination of human IL-2 in serum, plasma, and cell culture supernatant.

## Principle of the Test

The BD OptEIA test is a solid phase sandwich ELISA (Enzyme-Linked Immunosorbent Assay). It utilizes a monoclonal antibody specific for IL-2 coated on a 96-well plate. Standards and samples are added to the wells, and any IL-2 present binds to the immobilized antibody. The wells are washed and streptavidin-horseradish peroxidase conjugate mixed with biotinylated anti-human IL-2 antibody is added, producing an antibody-antigen-antibody “sandwich”. The wells are again washed and TMB substrate solution is added, which produces a blue color in direct proportion to the amount of IL-2 present in the initial sample. The Stop Solution changes the color from blue to yellow, and the microwell absorbances are read at 450 nm.

## Reagents Provided

Antibody Coated Wells:	2 plates of 96 breakable wells (12 strips × 8 wells) coated with anti-human IL-2 monoclonal antibody
Detection Antibody:	30 mL of biotinylated anti-human IL-2 monoclonal antibody with 0.15% ProClin®-150 as preservative
Standards:	4 vials lyophilized recombinant human IL-2
Enzyme Concentrate (250×):	150 µL of 250× concentrated Streptavidin-horseradish peroxidase conjugate with BSA* and ProClin®-150 as preservative
Standard/Sample Diluent:	30 mL of animal serum* with 0.09% sodium azide as preservative
ELISA Diluent:	12 mL of a buffered protein base with 0.09% sodium azide as preservative
Wash Concentrate (20×):	100 mL of 20× concentrated detergent solution with ProClin®-150 as preservative
TMB One-Step Substrate Reagent:	30 mL of 3,3',5,5'-tetramethylbenzidine (TMB) in buffered solution
Stop Solution:	13 mL of 1 M phosphoric acid
Plate Sealers:	4 sheets with adhesive backing

*\*Source of all serum proteins is from USDA inspected abattoirs located in the United States*

## Materials Required but not Provided

- Microplate reader capable of measuring absorbance at 450 nm
- Precision pipettes to deliver 50 µL and 100 µL volumes
- Adjustable 1 mL, 5 mL, 10 mL, 25 mL pipettes for reagent preparation
- Deionized or distilled water
- Wash bottle or automated microplate washer
- Log-log graph paper or computer and software for ELISA data analysis
- Tubes to prepare standard dilutions
- Laboratory timer
- Absorbent paper

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

# Storage Information

1. Store kit at 2 - 8°C. Do not use kit after expiration date.
2. Before use, bring all reagents to room temperature (18 - 25°C). Immediately after use, return to proper storage conditions.
3. Lyophilized standards are stable until kit expiration date. After reconstitution, use freshly reconstituted standard within 12 hours (stored at 2 - 8°C).

# Warnings and Precautions

1. Reagents that contain preservatives may be toxic if they are ingested, inhaled, or come in contact with skin.
2. Avoid contact of skin, eyes, or clothing with Stop Solution or Substrate Reagents.
3. Handle all serum and plasma specimens in accordance with NCCLS guidelines for preventing transmission of blood-borne infections.
4. Standard/Sample Diluent and ELISA Diluent contain less than 0.1% sodium azide. Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
5. **Warning**

Wash Concentrate (20X) (component 51-9003738) contains 0.002% (w/w), Human IL-2 Lyophilized Standard (component 51-27236E) contains 0.03% (w/w) and Detection Antibody Biotin Anti-Human IL-2 (component 51-27232E) contains 0.003% (w/w) of a CMIT/MIT mixture (3:1), which is a mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC No 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC No 220-239-6] (3:1).

## ***Hazard statements***

*May cause an allergic skin reaction.*

## ***Precautionary statements***

*Wear protective gloves / eye protection.*

*Wear protective clothing.*

*Avoid breathing mist/vapours/spray.*

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

*If skin irritation or rash occurs: Get medical advice/attention.*

*IF ON SKIN: Wash with plenty of water.*

*Dispose of contents/container in accordance with local/regional/national/international regulations.*

## **6. Danger**

Stop Solution (component 51-2608KC) contains 15.23% phosphoric acid (w/w).

### ***Hazard statements***

*Under normal conditions of intended use, this material is not expected to be an inhalation hazard.*

*Ingestion may cause severe irritation of the mouth, the esophagus and the gastrointestinal tract.*

*Causes severe skin burns and eye damage.*

### ***Precautionary statements***

*Wear protective gloves / eye protection.*

*Wear protective clothing.*

*IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.*

*IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.*

*Continue rinsing.*

*IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.*

*Dispose of contents/container in accordance with local/regional/national/international regulations.*

7. Enzyme Reagent (Streptavidin-horseradish peroxidase conjugate (SAv-HRP)) contains the following preservatives, 1,2-benzisothiazol-3(2H)-one (CAS: 2634-33-5) and CMIT/MIT 3:1 (CAS: 55965-84-9), which is a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC No 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC No 220-239-6] (3:1). The concentration of this material in the enzyme reagent is not classified as hazardous per GHS and CLP. Good chemical hygiene practices are recommended when handling any laboratory chemicals.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

# Specimen Collection and Handling

Specimens should be clear, non-hemolyzed and non-lipemic. Samples with expected values higher than the top standard, 500 pg/mL, should be diluted with Standard/Sample Diluent prior to running the assay.

**Cell culture supernatants:** Remove any particulate material by centrifugation and assay immediately or store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

**Serum:** Use a serum tube (eg, BD Vacutainer® Cat. No. 366430) and allow samples to clot for 30 minutes, then centrifuge for 10 minutes at 1000 x g. Remove serum and assay immediately or store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.


**Plasma:** Collect plasma using citrate, EDTA, or heparin as anticoagulant. Centrifuge for 10 minutes at 1000 x g within 30 minutes of collection. Assay immediately or store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

## Reagent Preparation

1. Bring all reagents to room temperature (18 - 25°C) before use.
2. Standards
  - a. Reconstitute 1 vial lyophilized Standard with required volume (noted on vial label) of Standard/Sample Diluent to prepare a 500 pg/mL stock standard. Allow the standard to equilibrate for at least 15 minutes before making dilutions. Gently vortex to mix.
  - b. Add 300  $\mu\text{L}$  Standard/Sample Diluent to 6 tubes. Label as 250 pg/mL, 125 pg/mL, 62.5 pg/mL, 31.3 pg/mL, 15.6 pg/mL, and 7.8 pg/mL.
  - c. Perform serial dilutions by adding 300  $\mu\text{L}$  of each standard to each standard to the next tube and vortexing between each transfer. The undiluted standard serves as the high standard



(500 pg/mL). The Standard/Sample Diluent serves as the zero standard (0 pg/mL).



### 3. Working Detector

**Note:** One-step incubation of Biotin/Avidin reagents.  
See *Assay Procedure, step 5*.

### 4. Wash Buffer

**Note:** If the Wash Concentrate contains visible crystals, warm to room temperature and mix gently until dissolved. Dilute required quantity of 20× Wash Concentrate with deionized or distilled water, mix. (To prepare 2,000 mL, add 100 mL Wash Concentrate to 1,900 mL water. At least 500 mL solution should be prepared for a full 96-well plate).

### 5. TMB One-Step Substrate Reagent

No more than 15 minutes prior to use, add required volume of TMB One-Step Substrate Reagent to a clean tube or reservoir. To prevent contamination, pipette out from the tube/reservoir instead of directly from bottle. Avoid prolonged exposure to light or contact with metal, air, or extreme temperature as color may develop.

# Assay Procedure

1. Bring all reagents and samples to room temperature (18 - 25°C) prior to use. It is recommended that all standards and samples be run in duplicate. A standard curve is required in each assay run.
2. Remove required quantity of test strips/wells, place in well holder.  
**Note:** Wells are provided in breakable 8-well strips. Strips may be “broken” into individual wells, replaced in well holder, and assayed. Return any unused wells to sealed pouch for 2 - 8°C storage.
3. Pipette 50 µL of ELISA Diluent into each well.
4. Pipette 100 µL of each standard (*see Reagent Preparation, step 2*) and sample into appropriate wells. Gently shake/tap the plate for 5 seconds to mix. Cover wells with Plate Sealer and incubate for 2 hours at room temperature.
5. Prepare Working Detector. Within 15 minutes prior to use, pipette required volume of Detection Antibody into a clean tube or flask. Add in required quantity of Enzyme Concentrate (250×), vortex or mix well. For a full 96-well plate, add 48 µL of Enzyme Concentrate into 12 mL of Detection Antibody.
6. Decant or aspirate contents of wells. Wash wells by filling with at least 300 µL/well prepared Wash Buffer (*see Reagent Preparation, step 4*), followed by decanting/aspirating. Repeat wash 4 times for a total of 5 washes. After the last wash, blot plate on absorbent paper to remove any residual buffer. Complete removal of liquid is required for proper performance.
7. Add 100 µL of prepared Working Detector (*see step 5 above*) to each well. Cover wells with Plate Sealer and incubate for 1 hour at room temperature.
8. Wash wells as in Step 6, but a total of 7 times.  
**Note:** In this final wash step, soak wells in wash buffer for 30 seconds to 1 minute for each wash. Thorough washing at this step is very important.
9. Add 100 µL of TMB One-Step Substrate Reagent to each well. Incubate plate (without Plate Sealer) for 30 minutes at room temperature in the dark.
10. Add 50 µL of Stop Solution to each well.
11. Read absorbance at 450 nm within 30 minutes of stopping reaction. If wavelength correction is available, subtract the optical density readings at 570 nm from readings at 450 nm.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

## Assay Procedure Summary

1. Add 50  $\mu$ L ELISA Diluent to each well.
2. Add 100  $\mu$ L standard or sample to each well.  
Incubate 2 hours at room temperature.
3. Aspirate and wash 5 times.
4. Add 100  $\mu$ L prepared Working Detector to each well.  
Incubate 1 hour at room temperature.
5. Aspirate and wash/soak 7 times.
6. Add 100  $\mu$ L TMB One-Step Substrate Reagent to each well.  
Incubate 30 minutes at room temperature.
7. Add 50  $\mu$ L Stop Solution to each well.  
Read at 450 nm within 30 minutes.  
 $\lambda$  correction 570 nm.

## Calculation of Results

Calculate the mean absorbance for each set of duplicate standards, controls and samples. Subtract the mean zero standard absorbance (ie, plate background) from each.


Plot the standard curve on log-log graph paper, with IL-2 concentration on the x-axis and absorbance on the y-axis. Draw the best fit straight line through the standard points.

To determine the IL-2 concentration of the unknowns, find the unknowns' mean absorbance value on the y-axis and draw a horizontal line to the standard curve. At the point of intersection, draw a vertical line to the x-axis and read the IL-2 concentration. If samples were diluted, multiply the interpolated IL-2 concentration by the dilution factor.

Computer-based curve-fitting statistical software may also be employed.

# Typical Data

This standard curve is for demonstration only. A standard curve must be run with each assay.



Concentration (pg/mL)	OD1	OD2	Mean	Zero Standard Subtracted
0	0.028	0.030	0.029	0.000
7.8	0.083	0.086	0.085	0.055
15.6	0.144	0.146	0.145	0.116
31.3	0.273	0.266	0.270	0.240
62.5	0.515	0.507	0.511	0.482
125	0.997	0.898	0.948	0.918
250	1.856	1.812	1.834	1.805
500	3.307	3.367	3.337	3.308

## Limitations of the Procedure

1. This kit is intended for use as an integral unit. Do not mix reagents from different kit lots. Reagents from other manufacturers/other available antibody clones should not be used in this kit.
2. Interference by drug metabolites, soluble receptors, or other binding proteins in specimens has not been thoroughly investigated. The possibility of interference cannot be excluded.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

# Performance

## Limit of Detection

The minimum detectable dose of IL-2 was determined to be 1 pg/mL. This is defined as two standard deviations above the mean optical density of 20 replicates of the zero standard.

## Recovery

Three different levels of IL-2 were spiked into samples of various matrices. Results are compared with the same amounts of IL-2 spiked into Standard/Sample Diluent, as follows:

	<b>Spike Concentration (pg/mL)</b>	<b>Average % Recovery</b>	<b>Range</b>
<b>Serum (n = 9)</b>	250	88	84 - 92
	125	86	70 - 107
	62.5	95	85 - 110
<b>Plasma (n = 5)</b>	250	83	77 - 91
	125	82	73 - 91
	62.5	90	78 - 100
<b>Cell culture media (n = 3)</b>	250	105	101 - 108
	125	116	110 - 120
	62.5	119	108 - 127

## Linearity

Samples spiked with high concentrations of IL-2 were serially diluted with Standard/Sample Diluent and run in the BD OptEIA Human IL-2 ELISA Kit II. Results are as follows:

<b>Dilution</b>		<b>Serum (n = 9)</b>	<b>Plasma (n = 5)</b>	<b>Cell culture media (n = 3)</b>
1:2	Average% of Expected Range	106 101 - 111	110 105 - 116	95 90 - 100
1:4	Average% of Expected Range	107 101 - 111	108 102 - 122	96 93 - 98
1:8	Average% of Expected Range	100 93 - 106	108 101 - 119	95 92 - 101
1:16	Average% of Expected Range	100 94 - 111	109 100 - 128	91 87 - 94

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

## Specificity

**Cross-Reactivity:** The factors listed below were spiked in Standard/Sample Diluent at 50 ng/mL to test for any cross-reactivity with the BD OptEIA™ Human IL-2 assay. No cross reactivity was identified.

### Recombinant Human

sCD23, Eotaxin, Fas, GM-CSF, GRO- $\alpha$ , GRO- $\beta$ , GRO- $\gamma$ , I-309, IFN- $\gamma$ , IL-3, IL-4, IL-5, IL-6, IL-8, IL-10, IL-12 (p40), IL-12(p70), IL-15, IL-16, IP-10, MCP-1, MCP-2, MCP-3, MIG, MIP-1 $\alpha$ , MIP-1 $\beta$ , NAP-2, PF-4, SDF-1 $\alpha$ , TNF- $\alpha$ , TNF- $\beta$

### Recombinant Mouse

IFN- $\gamma$ , IL-2, IL-3, IL-4, IL-5, IL-6, IL-9, IL-10, IL-12 (p40), IL-12 (p70), GM-CSF, MCP-1, MIG, MIP-1 $\alpha$ , MIP-1 $\beta$ , TCA-3, TNF- $\alpha$

### Recombinant Rat

IL-2, IL-10, GM-CSF, MCP-1, RANTES

**Interfering Substances:** The following substances at levels > 10 mg/mL were added to Standard/Sample Diluent spiked with 100 pg/mL IL-2. No effect on assay results was observed.

- Bilirubin
- Human hemoglobin
- Human transferrin
- Triglycerides
- Heparin (600 units/mL)
- Sodium Citrate
- EDTA

## Precision

### Intra-assay

Twenty-four replicates each of three different levels of IL-2 were tested in one plate. The following results were observed:

<b>Number of replicates</b>	24	24	24
<b>Mean Concentration</b>	233.7 pg/mL	111.0 pg/mL	56.4 pg/mL
<b>SD</b>	7.3	6.3	2.8
<b>%CV</b>	3.1	5.7	5.0

### Inter-assay

Three different levels of IL-2 were tested in four different plates. The following results were observed:

<b>Number of replicates</b>	32	32	32
<b>Mean Concentration</b>	238.2 pg/ml	119.9 pg/mL	59.6 pg/mL
<b>SD</b>	9.8	4.0	2.7
<b>%CV</b>	4.1	3.3	4.5

## Standardization

This immunoassay is calibrated against recombinant human IL-2.

The NIBSC/WHO First International Standard 86/504 (natural human IL-2) was evaluated in this kit. The conversion factor for NIBSC material is as follows:

NIBSC (86/504) equivalent value (IU/mL) = 0.0096 x BD OptEIA IL-2 value (pg/mL)

∴ 1 µg NIBSC IL-2 = 1.37 µg BD OptEIA recombinant human IL-2

## Experimental Results

### Serum

Twenty-six serum samples were tested in this assay. Nineteen samples measured less than the lowest standard, 7.8 pg/mL. Seven samples measured between 8.1 - 20.8 pg/mL.

## Plasma

Eight EDTA plasma samples were tested in this assay. All eight samples measured less than the lowest standard, 7.8 pg/mL.

## Cell Culture Supernatants

Human peripheral blood mononuclear cells or CD8<sup>+</sup> cells from five donors were stimulated with immobilized anti-human CD3, soluble anti-human CD28, recombinant human IL-2, and recombinant human IL-4 for 2 days. Cells were harvested and recultured in complete tissue culture medium containing rhIL-2 and rhIL-4 for 3 days. Finally, the *in-vitro* primed cells were harvested, washed, and restimulated with PMA and ionomycin for 4 hours, 24 hours, or 32 hours (as noted). Culture supernatants were collected and quantified for IL-2 using a BD OptEIA human IL-2 Kit II.

The results are as follows:

Donor No.	IL-2 (pg/mL)	
1	7,268.6	PBMC, 4 hr PMA/ionomycin stimulation
2	8,658.8	PBMC, 4 hr PMA/ionomycin stimulation
3	10,598.1	PBMC, 4 hr PMA/ionomycin stimulation
4	20,757.3	PBMC, 24 hr PMA/ionomycin stimulation
5	23,963.3	PBMC, 24 hr PMA/ionomycin stimulation
5	90,070.8	CD8 <sup>+</sup> cells, 32 hr PMA/ionomycin stimulation



# Troubleshooting

Problem	Possible Source	Corrective Action
Poor Precision	<ul style="list-style-type: none"> <li>• Inadequate washing / aspiration of wells</li> <li>• Inadequate mixing of reagents</li> <li>• Imprecise / inaccurate pipetting</li> <li>• Imprecise sealing of plate</li> </ul>	<ul style="list-style-type: none"> <li>• Check function of washing system</li> <li>• Ensure adequate mixing</li> <li>• Check / calibrate pipettes</li> <li>• Ensure complete sealing of plate</li> </ul>
Poor Standard Curve	<ul style="list-style-type: none"> <li>• Improper standard handling / dilution</li> <li>• Incomplete washing / aspiration of wells</li> <li>• Imprecise / inaccurate pipetting</li> </ul>	<ul style="list-style-type: none"> <li>• Ensure correct preparation of standards</li> <li>• Check function of washing system</li> <li>• Check / calibrate pipettes</li> </ul>
Low Signal	<ul style="list-style-type: none"> <li>• Inadequate reagent volumes added to wells</li> <li>• Incorrect incubation times / temperature</li> <li>• Overly high wash / aspiration pressure from automated plate-washer.</li> </ul>	<ul style="list-style-type: none"> <li>• Check / calibrate pipettes</li> <li>• Ensure sufficient incubation times / reagents warmed to room temperature</li> <li>• Utilize manual washing</li> </ul>

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

## References

1. Taniguchi, T., et al. 1983. Structure and expression of a cloned cDNA for human interleukin-2. *Nature*. 302:305-307.
2. Smith, K.A. 1988. Interleukin-2: inception, impact, and implications. *Science* 240:1169-1176.
3. Stern, A.S., et al. 1984. Purification to homogeneity and partial characterization of interleukin-2 from a human T-cell leukemia. *Proc. Natl. Acad. Sci. USA* 81:871-875.
4. Gillis, S., et al. 1978. T cell growth factor: Parameters of production and a quantitative microassay for activity. *J. Immunol.* 120:2027-2032.
5. Prussin, C., and D. D. Metcalfe. 1995. Detection of intracytoplasmic cytokine using flow cytometry and directly conjugated anti-cytokine antibodies. *J. Immunol. Meth.* 188:117-28.

# Plate Templates

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

[illegible]

## Notes

[illegible]

**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

## Notes

[illegible]

For Research Use Only. Not for use in diagnostic or therapeutic procedures.



**United States**

877.232.8995

**Canada**

866.979.9408

**Europe**

32.2.400.98.95

**Japan**

0120.8555.90

**Asia/Pacific**

65.66642770

**Latin America/Caribbean**

55.11.5185.9625

Toll Free 0800.771.71.57

**Becton, Dickinson and Company**

**BD Life Sciences - Biosciences**

2350 Qume Drive, San Jose, CA 95131

Ordering (US) 855.236.2772

Technical Service 877.232.8995

Fax 800.325.9637

[answers@bd.com](mailto:answers@bd.com)

[bdbiosciences.com](http://bdbiosciences.com)

Catalog No. 550611

Rev# 8

