

**BD OptEIA™**

# Human IFN- $\gamma$ ELISA Kit II

## Instruction Manual

Catalog No. 550612



*BD flow cytometers are Class 1 Laser Products.*

*For Research Use Only. Not for use in diagnostic or therapeutic procedures.*

*© 2018, Becton, Dickinson and Company. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means: electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without prior written permission from BD Biosciences.*

*The information in this guide is subject to change without notice. BD Biosciences reserves the right to change its products and services at any time to incorporate the latest technological developments. Although this guide has been prepared with every precaution to ensure accuracy, BD Biosciences assumes no liability for any errors or omissions, nor for any damages resulting from the application or use of this information.*

*BD Biosciences welcomes customer input on corrections and suggestions for improvement.*

*© 2018 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.*

# Table of Contents

Introduction.....	4
Principle of the Test.....	4
Reagents Provided.....	5
Materials Required but not Provided.....	5
Storage Information.....	6
Warnings and Precautions.....	6
Specimen Collection and Handling.....	8
Reagent Preparation.....	8
Assay Procedure.....	9
Assay Procedure Summary.....	10
Calculation of Results.....	11
Typical Data.....	11
Limitations of the Procedure.....	12
Performance.....	12
Limit of Detection.....	12
Recovery.....	13
Linearity.....	13
Specificity.....	13
Precision.....	15
Intra-assay.....	15
Inter-assay.....	15
Standardization.....	15
Experimental Results.....	16
Serum.....	16
Cell Culture Supernatants.....	16
Troubleshooting.....	17
References.....	18
Plate Templates.....	19

# Introduction

Interferon- $\gamma$  (IFN- $\gamma$ ) or immune interferon, is a potent multifunctional Type II interferon which is secreted by activated natural killer (NK) cells and T cells. IFN- $\gamma$  binds to a different receptor than the Type I viral interferons IFN- $\alpha$  and IFN- $\beta$ , with which little sequence homology is shared.

In addition to its antiviral effects, IFN- $\gamma$  can upregulate a number of lymphoid cell functions, including the antimicrobial and antitumor responses of macrophages, NK cells, and neutrophils. In addition, IFN- $\gamma$  can exert strong regulatory influences on the proliferation, differentiation, and effector responses of B cell and T cell subsets. These influences can involve IFN- $\gamma$ 's capacity to boost MHC class I and II expression by antigen-presenting cells as well as direct effects on B cells and T cells themselves. Human IFN- $\gamma$  is a 35 kD protein containing 143 amino acid residues. This cytokine also upregulates ICAM-1 expression on endothelial cells, and has become a primary marker for Th1-like immune response research.

The BD OptEIA™ ELISA Kit II format was developed for superior accuracy with serum and plasma specimens. The data that demonstrates this enhancement can be located in the Performance “Recovery” and “Linearity” sections.

The BD OptEIA Human IFN- $\gamma$  ELISA Kit II is for the *in vitro* quantitative determination of human IFN- $\gamma$  in serum, plasma, and cell culture supernatant.

## Principle of the Test

The BD OptEIA test is a solid phase sandwich ELISA (Enzyme-Linked Immunosorbent Assay). It utilizes a monoclonal antibody specific for IFN- $\gamma$  coated on a 96-well plate. Standards and samples are added to the wells, and any IFN- $\gamma$  present binds to the immobilized antibody. The wells are washed and streptavidin-horseradish peroxidase conjugate mixed with biotinylated anti-human IFN- $\gamma$  antibody is added, producing an antibody-antigen-antibody “sandwich”. The wells are again washed and TMB substrate solution is added, which produces a blue color in direct proportion to the amount of IFN- $\gamma$  present in the initial sample. The Stop Solution changes the color from blue to yellow, and the microwell absorbances are read at 450 nm.

## Reagents Provided

Antibody Coated Wells:	2 plates of 96 breakable wells (12 strips × 8 wells) coated with anti-human IFN- $\gamma$ monoclonal antibody
Detection Antibody:	30 mL of biotinylated anti-human IFN- $\gamma$ monoclonal antibody with ProClin™- 150 as preservative
Standards:	4 vials lyophilized recombinant human IFN- $\gamma$
Enzyme Concentrate (250×):	150 $\mu$ L of 250× concentrated Streptavidin-horseradish peroxidase conjugate with BSA* and ProClin™- 150 as preservative
Standard/Sample Diluent:	30 mL of animal serum* with 0.09% sodium azide as preservative
ELISA Diluent:	12 mL of a buffered protein base with 0.09% sodium azide as preservative
Wash Concentrate (20×):	100 mL of 20× concentrated detergent solution with ProClin™- 150 as preservative
TMB One-Step Substrate Reagent:	30 mL of 3,3',5,5'-tetramethylbenzidine (TMB) in buffered solution
Stop Solution:	13 mL of 1 M phosphoric acid
Plate Sealers:	4 sheets with adhesive backing

*\*Source of all serum proteins is from USDA inspected abattoirs located in the United States*

## Materials Required but not Provided

- Microplate reader capable of measuring absorbance at 450 nm
- Precision pipettes to deliver 50  $\mu$ L and 100  $\mu$ L volumes
- Adjustable 1 mL, 5 mL, 10 mL, 25 mL pipettes for reagent preparation
- Deionized or distilled water
- Wash bottle or automated microplate washer
- Log-log graph paper or computer and software for ELISA data analysis
- Tubes to prepare standard dilutions
- Laboratory timer
- Absorbent paper

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

# Storage Information

1. Store kit at 2 - 8°C. Do not use kit after expiration date.
2. Before use, bring all reagents to room temperature (18 - 25°C). Immediately after use, return to proper storage conditions.
3. Lyophilized standards are stable until kit expiration date. After reconstitution, use freshly reconstituted standard within 12 hours (stored at 2 - 8°C).

# Warnings and Precautions

1. Reagents that contain preservatives may be toxic if ingested, inhaled, or brought in contact with skin.
2. Avoid contact of skin, eyes, or clothing with Stop Solution or Substrate Reagents.
3. Handle all serum and plasma specimens in accordance with NCCLS guidelines for preventing transmission of blood-borne infections.
4. Standard/Sample Diluent and ELISA Diluent contain less than 0.1% sodium azide. Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
5. **Warning**

Wash Concentrate (20X) (component 51-9003738) contains 0.002% (w/w), Human IFN-gamma Lyophilized Standard (component 51-27226E) contains 0.03% (w/w) and Detection Antibody Biotin Anti-Human IFN-gamma (component 51-27222E) contains 0.003% (w/w) of a CMIT/MIT mixture (3:1), which is a mixture of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC No 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC No 220-239-6] (3:1).

## **Hazard statements**

*May cause an allergic skin reaction.*

## **Precautionary statements**

*Wear protective gloves / eye protection.*

*Wear protective clothing.*

*Avoid breathing mist/vapours/spray.*

*If skin irritation or rash occurs: Get medical advice/attention.*

**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

*IF ON SKIN: Wash with plenty of water.*

*Dispose of contents/container in accordance with local/regional/national/international regulations.*

## **6. Danger**

Stop Solution (component 51-2608KC) contains 15.23% phosphoric acid (w/w).

### **Hazard statements**

*Under normal conditions of intended use, this material is not expected to be an inhalation hazard.*

*Ingestion may cause severe irritation of the mouth, the esophagus and the gastrointestinal tract.*

*Causes severe skin burns and eye damage.*

### **Precautionary statements**

*Wear protective gloves / eye protection.*

*Wear protective clothing.*

*IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.*

*IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.*

*Continue rinsing.*

*IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.*

*Dispose of contents/container in accordance with local/regional/national/international regulations.*

7. Enzyme Reagent (Streptavidin-horseradish peroxidase conjugate (SAv-HRP)) contains the following preservatives, 1,2-benzisothiazol-3(2H)-one (CAS: 2634-33-5) and CMIT/MIT 3:1 (CAS: 55965-84-9), which is a mixture of 5-chloro-2-methyl-4-isothiazolin-3-one [EC No 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC No 220-239-6] (3:1). The concentration of this material in the enzyme reagent is not classified as hazardous per GHS and CLP. Good chemical hygiene practices are recommended when handling any laboratory chemicals.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

# Specimen Collection and Handling

Specimens should be clear, non-hemolyzed and non-lipemic. Samples with expected values higher than the top standard, 300 pg/mL, should be diluted with Standard/Sample Diluent prior to running the assay.

**Cell culture supernatants:** Remove any particulate material by centrifugation and assay immediately or store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

**Serum:** Use a serum tube (eg, BD Vacutainer® Cat. No. 366430) and allow samples to clot for 30 minutes, then centrifuge for 10 minutes at  $1000 \times g$ . Remove serum and assay immediately or store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.


**Plasma:** Collect plasma using citrate, EDTA, or heparin as anticoagulant. Centrifuge for 10 minutes at  $1000 \times g$  within 30 minutes of collection. Assay immediately or store samples at  $\leq -20^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles.

## Reagent Preparation

1. Bring all reagents to room temperature ( $18 - 25^{\circ}\text{C}$ ) before use.
2. Standards
  - a. Reconstitute 1 vial lyophilized Standard with required volume (noted on vial label) of Standard/Sample Diluent to prepare a 300 pg/mL stock standard. Allow the standard to equilibrate for at least 15 minutes before making dilutions. Gently vortex to mix.
  - b. Add 300  $\mu\text{L}$  Standard/Sample Diluent to 6 tubes. Label as 150 pg/mL, 75 pg/mL, 37.5 pg/mL, 18.8 pg/mL, 9.4 pg/mL, and 4.7 pg/mL.
  - c. Perform serial dilutions by adding 300  $\mu\text{L}$  of each standard to the next tube and vortexing between each transfer. The undiluted



standard serves as the high standard (300 pg/mL). The Standard/ Sample Diluent serves as the zero standard (0 pg/mL).



### 3. Working Detector

**Note:** One-step incubation of Biotin/Streptavidin reagents.  
See *Assay Procedure*, step 5.

### 4. Wash Buffer

**Note:** If the Wash Concentrate contains visible crystals, warm to room temperature and mix gently until dissolved. Dilute required quantity of 20× Wash Concentrate with deionized or distilled water, mix. (To prepare 2,000 mL, add 100 mL Wash Concentrate to 1,900 mL water. At least 500 mL solution should be prepared for a full 96-well plate).

### 5. TMB One-Step Substrate Reagent

No more than 15 minutes prior to use, add required volume of TMB One-Step Substrate Reagent to a clean tube or reservoir. To prevent contamination, pipette out from the tube/ reservoir instead of directly from bottle. Avoid prolonged exposure to light or contact with metal, air, or extreme temperature as color may develop.

## Assay Procedure

1. Bring all reagents and samples to room temperature (18 - 25°C) prior to use. It is recommended that all standards and samples be run in duplicate. A standard curve is required in each assay run.

2. Remove required quantity of test strips/wells, place in well holder.

**Note:** Wells are provided in breakable 8-well strips. Strips may be “broken” into individual wells, replaced in well holder, and assayed. Return any unused wells to sealed pouch for 2 - 8°C storage.

3. Pipette 50 µL of ELISA Diluent into each well.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

4. Pipette 100  $\mu$ L of each standard (see *Reagent Preparation*, step 2) and sample into appropriate wells. Gently shake/tap the plate for 5 seconds to mix. Cover wells with Plate Sealer and incubate for 2 hours at room temperature.
5. Prepare Working Detector. Within 15 minutes prior to use, pipette required volume of Detection Antibody into a clean tube or flask. Add in required quantity of Enzyme Concentrate (250 $\times$ ), vortex or mix well. For a full 96-well plate, add 48  $\mu$ L of Enzyme Concentrate into 12 mL of Detection Antibody.
6. Decant or aspirate contents of wells. Wash wells by filling with at least 300  $\mu$ L/well prepared Wash Buffer (see *Reagent Preparation*, step 4), followed by decanting/aspirating. Repeat wash 4 times for a total of 5 washes. After the last wash, blot plate on absorbent paper to remove any residual buffer. Complete removal of liquid is required for proper performance.
7. Add 100  $\mu$ L of prepared Working Detector (see *step 5* above) to each well. Cover wells with Plate Sealer and incubate for 1 hour at room temperature.
8. Wash wells as in Step 6, but a total of 7 times.  
**Note:** In this final wash step, soak wells in wash buffer for 30 seconds to 1 minute for each wash. Thorough washing at this step is very important.
9. Add 100  $\mu$ L of TMB One-Step Substrate Reagent to each well. Incubate plate (without Plate Sealer) for 30 minutes at room temperature in the dark.
10. Add 50  $\mu$ L of Stop Solution to each well.
11. Read absorbance at 450 nm within 30 minutes of stopping reaction. If wavelength correction is available, subtract the optical density readings at 570 nm from readings at 450 nm.

## Assay Procedure Summary

1. Add 50  $\mu$ L ELISA Diluent to each well.
2. Add 100  $\mu$ L standard or sample to each well.  
Incubate 2 hours at room temperature.
3. Aspirate and wash 5 times.
4. Add 100  $\mu$ L prepared Working Detector to each well.  
Incubate 1 hour at room temperature.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

5. Aspirate and wash/soak 7 times.
6. Add 100  $\mu$ L TMB One-Step Substrate Reagent to each well.  
Incubate 30 minutes at room temperature.
7. Add 50  $\mu$ L Stop Solution to each well.  
Read at 450 nm within 30 minutes.  
 $\lambda$  correction 570 nm.

## Calculation of Results

Calculate the mean absorbance for each set of duplicate standards, controls and samples. Subtract the mean zero standard absorbance (ie, plate background) from each.


Plot the standard curve on log-log graph paper, with IFN- $\gamma$  concentration on the x-axis and absorbance on the y-axis. Draw the best fit straight line through the standard points.

To determine the IFN- $\gamma$  concentration of the unknowns, find the unknowns' mean absorbance value on the y-axis and draw a horizontal line to the standard curve. At the point of intersection, draw a vertical line to the x-axis and read the IFN- $\gamma$  concentration. If samples were diluted, multiply the interpolated IFN- $\gamma$  concentration by the dilution factor.

Computer-based curve-fitting statistical software may also be employed.

## Typical Data

This standard curve is for demonstration only. A standard curve must be run with each assay.



For Research Use Only. Not for use in diagnostic or therapeutic procedures.

<b>Concentration (pg/mL)</b>	<b>OD1</b>	<b>OD2</b>	<b>Mean</b>	<b>Zero Standard Subtracted</b>
0	0.037	0.033	0.035	0.000
4.7	0.070	0.069	0.070	0.035
9.4	0.116	0.112	0.114	0.079
18.8	0.191	0.180	0.186	0.151
37.5	0.348	0.327	0.338	0.303
75	0.604	0.604	0.604	0.569
150	1.188	1.159	1.174	1.139
300	2.402	2.402	2.402	2.367

## Limitations of the Procedure

1. This kit is intended for use as an integral unit. Do not mix reagents from different kit lots. Reagents from other manufacturers/other available antibody clones should not be used in this kit.
2. Interference by drug metabolites, soluble receptors, or other binding proteins in specimens has not been thoroughly investigated. The possibility of interference cannot be excluded.

## Performance

### Limit of Detection

The minimum detectable dose of IFN- $\gamma$  was determined to be 1.0 pg/mL. This is defined as two standard deviations above the mean optical density of 20 replicates of the zero standard.

## Recovery

Three different levels of IFN- $\gamma$  were spiked into samples of various matrices, as noted below. Results are compared with the same amounts of IFN- $\gamma$  spiked into Standard/Sample Diluent, as follows:

	<b>Spike Concentration (pg/mL)</b>	<b>Average % Recovery</b>	<b>Range</b>
Serum (n = 9)	150	97	89 - 102
	75	97	89 - 110
	37.5	106	81 - 139
Plasma (n = 9)	150	94	88 - 101
	75	90	82 - 94
	37.5	89	76 - 97
Cell culture media (n = 3)	150	105	103 - 107
	75	105	100 - 114
	37.5	105	98 - 111

## Linearity

Various samples spiked with high concentrations of IFN- $\gamma$  were serially diluted with Standard/Sample Diluent and run in the BD OptEIA Human IFN- $\gamma$  ELISA Kit II. Results are as follows:

<b>Dilution</b>		<b>Serum (n = 9)</b>	<b>Plasma (n = 9)</b>	<b>Cell culture media (n = 3)</b>
1:2	Average % of Expected Range	96 91 - 105	100 96 - 103	98 94- 98
1:4	Average % of Expected Range	98 87 - 109	100 95 - 103	98 94 - 101
1:8	Average % of Expected Range	99 91 - 124	99 91 - 106	98 94 - 102
1:16	Average % of Expected Range	102 90 - 135	94 86 - 104	96 88 - 101

## Specificity

**Cross-Reactivity:** The proteins listed below were spiked in Standard/Sample Diluent at 100 ng/mL to test for any cross-reactivity with the BD OptEIA™ Human IFN- $\gamma$  assay. No cross reactivity was identified.

**Interference:** The factors listed below were spiked at 100 ng/mL in Standard/Sample Diluent with 100 pg/mL IFN- $\gamma$  to test for any interference with the quantitation of human IFN- $\gamma$ . No effect on assay results was observed.

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

## **Recombinant Human**

IL-1 $\alpha$ , IL-1 $\beta$ , IL-2, IL-3, IL-4, IL-5, IL-6, IL-7, IL-8, IL-9, IL-10, IL-11, IL-12 (p40), IL-12 (p70), IL-13, IL-15, TNF, Eotaxin, G-CSF (10 ng/mL), GM-CSF, GRO, CD23, Lymphotoxin (10 ng/mL), MIP-1 $\beta$ , MCP-1, MCP-2, MCP-3, MCP-4, NAP2, IP-10, NT-3, PDGF-AA, SCF (10 ng/mL), LT- $\alpha$ , VEGF

## **Recombinant Mouse**

IL-1 $\beta$ , IL-2, IL-3, IL-4, IL-5, IL-6, IL-7, IL-9, IL-10, IL-12 (p70), IL-15, IFN- $\gamma$ , GM-CSF, MCP-1, MIG, TCA3, TNF

## **Recombinant Rat**

IL-2, IL-4, IL-6, IL-10, GM-CSF, IFN- $\gamma$ , IL-12 (p70), RANTES

## **Other:**

Viral IL-10 (10 ng/mL)

## **Interfering Substances:**

The following substances at levels > 20 mg/mL were added to Standard/ Sample Diluent spiked with 200 pg/mL IFN- $\gamma$ . No effect on assay results was observed.

- Bilirubin
- Human hemoglobin
- Human transferrin
- Triglycerides
- Heparin (300 units/mL)
- Sodium Citrate
- EDTA

# Precision

## Intra-assay

Twenty-four replicates each of three different levels of IFN- $\gamma$  were tested in one plate. The following results were observed:

Number of Replicates	24	24	24
Mean Concentration	126.7 pg/mL	64.2 pg/mL	33.0 pg/mL
SD	4.0	2.7	1.3
%CV	3.2	4.1	4.0

## Inter-assay

Three different levels of IFN- $\gamma$  were tested in four different plates. The following results were observed:

Number of Replicates	32	32	32
Mean Concentration	132.3 pg/mL	68.1 pg/mL	35.7 pg/mL
SD	5.7	4.0	3.5
%CV	4.3	5.8	9.9

# Standardization

The Human IFN- $\gamma$  immunoassay is calibrated against recombinant human IFN- $\gamma$ .

The NIBSC/WHO First British Standard 82/587 (recombinant human IFN- $\gamma$  prepared by induction with phytohaemagglutinin) was evaluated in this kit. The conversion factor for NIBSC material is as follows:

NIBSC equivalent value (IU/mL) = 0.010  $\times$  BD OptEIA IFN- $\gamma$  value (pg/mL)

The NIAID Standard #Gxg01-902-535 (recombinant human IFN- $\gamma$  expressed in E. Coli) was evaluated in this kit. The conversion factor for NIAID material is as follows:

NIAID equivalent value (IU/mL) = 0.021  $\times$  BD OptEIA IFN- $\gamma$  value (pg/mL)

# Experimental Results

## Serum

Twenty-one serum samples were tested in the human IFN- $\gamma$  assay. Eighteen samples measured less than 4.7 pg/mL (lowest standard level). Three samples measured between 5.5 – 10.3 pg/mL.

## Cell Culture Supernatants

Human peripheral blood mononuclear cells or CD8<sup>+</sup> cells were stimulated with immobilized anti-human CD3, soluble anti-human CD28, recombinant human IL-2, and recombinant human IL-4 for 2 days. Cells were harvested and recultured in complete tissue culture medium containing rhIL-2 and rhIL-4 for 3 days. Finally, the *in vitro* primed cells were harvested, washed, and restimulated with PMA and ionomycin for 4 hours, 24 hours, or 32 hours (as noted). Culture supernatants were collected and quantified for IFN- $\gamma$  using a BD OptEIA Human IFN- $\gamma$  Kit II. The results are as follows:

Donor No.	IFN- $\gamma$ (pg/mL)	
1	703.9	PBMC, 4 hr PMA/ionomycin stimulation
2	4,556.2	PBMC, 4 hr PMA/ionomycin stimulation
3	2,665.6	PBMC, 4 hr PMA/ionomycin stimulation
4	4,490.6	PBMC, 24 hr PMA/ionomycin stimulation
5	5,043.6	PBMC, 24 hr PMA/ionomycin stimulation
5	88,941.0	CD8 <sup>+</sup> cells, 32 hr PMA/ionomycin stimulation



# Troubleshooting

Problem	Possible Source	Corrective Action
Poor Precision	<ul style="list-style-type: none"> <li>• Inadequate washing / aspiration of wells</li> <li>• Inadequate mixing of reagents</li> <li>• Imprecise / inaccurate pipetting</li> <li>• Imprecise sealing of plate</li> </ul>	<ul style="list-style-type: none"> <li>• Check function of washing system</li> <li>• Ensure adequate mixing</li> <li>• Check / calibrate pipettes</li> <li>• Ensure complete sealing of plate</li> </ul>
Poor Standard Curve	<ul style="list-style-type: none"> <li>• Improper standard handling / dilution</li> <li>• Incomplete washing / aspiration of wells</li> <li>• Imprecise / inaccurate pipetting</li> </ul>	<ul style="list-style-type: none"> <li>• Ensure correct preparation of standards</li> <li>• Check function of washing system</li> <li>• Check / calibrate pipettes</li> </ul>
Low Signal	<ul style="list-style-type: none"> <li>• Inadequate reagent volumes added to wells</li> <li>• Incorrect incubation times / temperature</li> <li>• Overly high wash / aspiration pressure from automated plate-washer.</li> </ul>	<ul style="list-style-type: none"> <li>• Check / calibrate pipettes</li> <li>• Ensure sufficient incubation times / reagents warmed to room temperature</li> <li>• Utilize manual washing</li> </ul>

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

## References

1. Farrar, M.A. and R.D. Schreiber (1993) The biology and biochemistry of interferon-gamma and its receptor. *Annu. Rev. Immunol.* 11: 571 - 611.
2. Grossberg, S.E. *et al.* (1989) Interferon receptors and their role in interferon action. *Experientia* 45: 508.
3. Gray, P.W. *et al.* (1982) Structure of the human immune interferon gene. *Nature* 295: 503.
4. Rinderknecht, E. *et al.* (1984) Natural human interferon-gamma. Complete amino acid sequence and determination of sites of glycosylation *J. Biol. Chem.* 259: 6790.
5. DeGrado, W.F. *et al.* (1982) Sequence and structural homologies among type I and type II interferons. *Nature* 300: 379.
6. Zoon, K.C. (1987) Human interferons: structure and function. *Interferon* 9: 1.

# Plate Templates

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

For Research Use Only. Not for use in diagnostic or therapeutic procedures.

[illegible]

## Notes

[illegible]

**For Research Use Only. Not for use in diagnostic or therapeutic procedures.**

## Notes

[illegible]



**United States**

877.232.8995

**Canada**

866.979.9408

**Europe**

32.2.400.98.95

**Japan**

0120.8555.90

**Asia/Pacific**

65.66642770

**Latin America/Caribbean**

55.11.5185.9625

Toll Free 0800.771.71.57

**Becton, Dickinson and Company**  
**BD Life Sciences - Biosciences**  
2350 Qume Drive, San Jose, CA 95131  
Ordering (US) 855.236.2772  
Technical Service 877.232.8995  
Fax 800.325.9637  
[answers@bd.com](mailto:answers@bd.com)

[bdbiosciences.com](http://bdbiosciences.com)

Catalog No. 550612

Rev# 7

