


# BD Falcon™ Pipetter Assurance Program



## An unprecedented service offer available to BD Falcon™ Pipet customers...

BD Biosciences is pleased to provide this free maintenance program to help ensure that all BD Falcon™ Pipetters in your laboratory are drip-free and in excellent working condition.

Upon registration, you will automatically receive two new rubber inserts for the nosepiece housing of each BD Falcon Pipetter in your laboratory at the beginning of every calendar year. The rubber insert is responsible for forming the contact seal between the BD Falcon pipetter nosepiece housing and the pipet itself. These inserts can wear and become inflexible with use. Routine replacement ensures a secure seal. We recommend replacing them every six months.

To register, fax the following information to **BD Biosciences** - Discovery Labware Technical Service Department at **978.901.7491**.

Name \_\_\_\_\_

Institution \_\_\_\_\_

Department \_\_\_\_\_ Bldg/Room# \_\_\_\_\_

Address \_\_\_\_\_

City \_\_\_\_\_ State \_\_\_\_\_ Zip Code \_\_\_\_\_

Telephone \_\_\_\_\_ Fax \_\_\_\_\_

**Enter below the number of pipettors to be registered:**

Portable or Corded Mechanical Pipetter(s)