

Maintenance Log for BD FACSCanto™ II System

Serial No. _____

Month/Year: _____

Initials and Dates																												
DAILY																												
DATE																												
Fluidics Startup																												
QC																												
7-color setup beads																												
P/F																												
CS&T																												
P/F																												
Export Experiments																												
Empty Condensation Trap																												
Cytometer Shutdown 1. Decontaminate SIT with 10% bleach solution (eg, BD FACSClean) and rinse with DI water. 2. Run Fluidics Shutdown.	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	/	
WEEKLY																												
Purge fluidics filters																												
Clean flow cell																												
Clean tube sensor																												
Archive experiments and FCS files																												
Back up database																												
Defragment hard drive																												
MONTHLY																												
Perform Long Clean																												
Change waste container cap																												

Maintenance Log for BD FACSCanto™ II System

Serial No. _____

Month/Year: _____

Initials and Dates

EVERY 6 MONTHS						
Perform Periodic Maintenance (PM)						
Replace air filters (or as needed)						
Replace fluidics filters (or as needed)						
Define Baseline						
Check application settings if baseline was changed						

UNSCHEDULED MAINTENANCE (Replace bal seal, fuses, fluidic sensors, etc)	Comments

LOADER MAINTENANCE																									
Clean optical sensors on Loader drawer																									