

FACSComp™

Advanced software from Becton Dickinson for automated instrument setup and monitoring

Now Becton Dickinson offers the

complete clinical solution for your

immunophenotyping needs.

The comprehensive FACS MultiSET™


System integrates powerful

state-of-the-art software and reagents

with the FACSCalibur™ flow cytometer

to bring precision to your counts and

improved productivity to your lab.




Today's clinical laboratory environment is characterized by increasing workloads and demands for improving productivity, while meeting strict quality assurance requirements. To help your laboratory meet these needs, FACSComp provides the convenience of automated daily instrument setup and also helps you track instrument performance. FACSComp supports automated setup for a broad range of assays, now including three- and four-color¹ applications such as TriTEST™ and MultiTEST™¹. Designed for use with Becton Dickinson's FACSCalibur, FACSort™, and FACScan™ instruments with FACStation™ systems, the latest evolution of FACSComp is an integral component of our fully supported systems and applications for two-, three-, and four-color^{1,2} flow cytometry.




FACSComp™

Advanced software from Becton Dickinson for automated instrument setup and monitoring



Convenient FACSComp Summary Report summarizes instrument settings and performance evaluation results.



In addition to providing automated instrument setup for a broad range of applications, FACSComp now includes an Optimization feature, allowing use of actual cell samples to confirm that instrument settings are appropriate for your specific applications.



Mac™ OS

Meeting Your Requirements for Daily Quality Assurance

FACSComp, when used with CaliBRITE™ beads, provides convenient routine instrument setup for immunophenotyping applications, as well as monitoring of instrument performance to ensure reproducible results—day after day. Instrument settings and performance test results from multiple runs can be saved automatically to an export file which can then be read by third-party spreadsheet or statistics software to generate Levey-Jennings plots.

Now Supports Four-Color Applications¹

The latest version of FACSComp software is a key component of the FACS MultiSET System, providing automated instrument setup for two-, three- and four-color applications. This feature will be essential for an increasing number of laboratories as new multicolor applications become widespread.

FACSComp offers setup options for both conventional lyse/wash methods (two-,³ three-,¹ and four-color¹) and Becton Dickinson's convenient lyse/no-wash methods (three-,³ and four-color¹) used with reagent systems such as MultiTEST,¹ TriTEST, FASTIMMUNE™,¹ and ProCOUNT™¹. By greatly reducing the time and effort required for your daily instrument setup, FACSComp allows you to focus on running samples and reporting results, not on the often cumbersome details of manual setup procedures for multicolor assays. FACSComp allows your laboratory to take full advantage of the productivity gains associated with MultiSET and other multicolor applications.

Leading the Way in Cell Analysis

For more information on how FACSComp can help you meet your needs for quality assurance and increased productivity, contact your local Becton Dickinson Immunocytometry Systems representative today.

¹ For research use only. Not for use in diagnostic or therapeutic procedures.

² FACSCalibur or FACSsort equipped with the 4-color option.

³ For in vitro diagnostic use.

©1997 Becton Dickinson and Company

FACSComp, FACStation, FACSCalibur, FACSsort, FACScan, MultiSET, TriTEST, FASTIMMUNE, MultiTEST, ProCOUNT, and CaliBRITE are trademarks of Becton Dickinson and Company. FACService is a service mark of Becton Dickinson and Company.

Macintosh is a registered trademark of Apple Computer Inc. Mac OS and Power Macintosh are trademarks of Apple Computer Inc.

**BECTON
DICKINSON**

Becton Dickinson
Immunocytometry Systems
2350 Qume Drive
San Jose, CA 95131-1807
Ordering Information
(800) 223-8226
Customer Support Center
(800) 448-BDIS (2347)
Fax (408) 954-BDIS (2347)

Becton Dickinson
Canada, Inc.
2464 South Sheridan Way
Mississauga, Ontario
L5J 2M8
Canada
Tel (905) 822-4820
Fax (905) 822-2644

Becton Dickinson
European HQ
Denderstraat 24
B-9320 Erembodegem-Aalst
Belgium
Tel (32) 53-720211
Fax (32) 53-720450

Nippon Becton Dickinson
Company, Ltd.
DS Bldg.
5-26, Akasaka 8-chome
Minato-ku, Tokyo 107
Japan
Tel (81) 3-5413-8251
Fax (81) 3-5413-8155

Becton Dickinson
Worldwide, Inc.
30 Tuas Avenue, #2
Singapore 2263
Tel (65) 861-0633
Fax (65) 860-1590

23-2858-01
1/97 10k



FM 32438