

Technical Data Sheet

Annexin V Binding Buffer, 10X concentrate

Product Information

Material Number: 556454

Size: 50 ml

Description

The Annexin V Binding Buffer, 10X concentrate has been formulated for use in flow cytometric assays utilizing Annexin V. The Annexin V Binding Buffer is a 10X concentrate composed of a 0.2 µm sterile filtered 0.1M Hepes (pH 7.4), 1.4M NaCl, and 25 mM CaCl₂ solution. Prior to staining cells, an appropriate quantity of a 1X working solution should be made by diluting the 10X concentrate 1:10 with distilled water. Any remaining 1X solution at the end of the experiment should be discarded.

Preparation and Storage

Store undiluted at 4° C.

Application Notes

Application

Flow cytometry	Routinely Tested
----------------	------------------

Suggested Companion Products

Catalog Number	Name	Size	Clone
556419	FITC Annexin V	200 tests	(none)
556421	PE Annexin V	200 tests	(none)
559934	Cy5 Annexin V	200 tests	(none)
556417	Biotin Annexin V	200 tests	(none)
550475	APC Annexin V	200 tests	(none)
559935	Cy5.5 Annexin V	100 tests	(none)

Product Notices

- Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

References

Vermes I, Haanen C, Steffens-Nakken H, Reutelingsperger C. A novel assay for apoptosis. Flow cytometric detection of phosphatidylserine expression on early apoptotic cells using fluorescein labelled Annexin V. *J Immunol Methods*. 1995; 184(1):39-51.(Methodology: Apoptosis, Flow cytometry)

BD Biosciences

bdbiosciences.com

United States 877.232.8995 Canada 888.259.0187 Europe 32.53.720.550 Japan 0120.8555.90 Asia Pacific 65.6861.0633 Latin America/Caribbean 55.11.5185.9995

For country-specific contact information, visit bdbiosciences.com/how_to_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2006 BD



BD

BD Biosciences