

BD FACSDuetTM sample preparation system

Product list



BD FACSDuetTM sample preparation system

Sample preparation system

Item	Description	Cat. No.
BD FACSDuet™ Sample Prep. System	Sample preparation system for all Lyse-No-Wash applications. Includes kit for physical integration with the BD FACSLyric™ flow cytometer. Also includes a starter accessory kit:	662588
	Reagent Loading Rack (Set of 2 racks)	
	Specimen Loading Rack (Set of 4 racks)	
	13x75 Specimen Tube Adapters (4 packs of 10 ea.)	
	BD FACSLyric Secondary Tube Rack (30 Tube Rack) (Set of 4 racks)	
	Handheld Barcode Scanner with Stand	
	Lysing Solution 1L Tank	
	10L DI Water Tank	
	10L Waste Tank	
	10L Saline Tank	
	50mL Bleach Container	
	50mL Chilled Saline Container	
	Documentation Kit (IFU, User Guide, SPG, SLG)	
	Panel PC	

Optional components

Item	Description	Cat. No.
Table (recommended when BD FACSDuet and BD FACSLyric are a physically integrated system)	200.0cm L x 97.0cm D x 86.35cm H Table includes BD FACSLyric flat screen desk support, dripping plate with receptacles, 2 chariots for BD FACSDuet and BD FACSLyric tanks	663375
Secondary Tube Labels Printer	Printer generating labels for 12x75mm secondary tubes	662609
Handheld Barcode Scanner, with Stand	External barcode reader for both BD FACSLyric and BD FACSDuet	663058
Uninterruptible Power Manager (UPM)	1440VA 230V IN/OUT	663407
Panel PC	Replacement touch screen PC	663211

BD FACSDuet[™] sample preparation system

Accessories

Item	Description	Cat. No.
Specimen Loading Rack		
Specimen Loading Rack (Set of 4)	Set of 4 racks for loading primary specimen tubes. Does not include tube adaptors	662592
Streck Cyto-Chex specimen Loading Rack (Set of 4)	Set of 4 racks for loading Streck Cyto-Chex tubes.	662593
Specimen Rack Adapters		
13x65 Specimen Tube Adapters (Set of 10)	13mmx65mm adapter for Sarstedt 2.6mL and 3.4mL tubes	662594
11x66 Specimen Tube Adapters (Set of 10)	11mmx66mm adapter for Sarstedt 2.7mL tubes	662595
13x75 Specimen Tube Adapters (Set of 10)	13mmx75mm adapter for BD 2mL, 3mL and 4mL, Greiner 2mL and 3mL, Streck 5mL tubes	662596
15x75 Specimen Tube Adapters (Set of 10)	15mmx75mm adapter for Sarstedt 4mL tubes	662597
13x90 Specimen Tube Adapters (Set of 10)	13mmx90mm adapter for Sarstedt 4.9mL tubes	662598
13x100 Specimen Tube Adapters (Set of 10)	13mmx100mm adapter for BD and Greiner 6mL and 7mL tubes	662600
16x100 Specimen Tube Adapters (Set of 10)	16mmx100mm adapter for BD 9.5mL and 10mL and Greiner 9mL and 10mL tubes	662599
Secondary Tube Racks		
BD FACSLyric Secondary Tube Rack (30 Tube Rack)	Secondary tube rack to accommodate 30 12x75mm tubes	651314
BD FACSLyric Secondary Tube Rack (40 Tube Rack)	Secondary tube rack to accommodate 40 12x75mm tubes	651319
Reagent Racks		
Reagent Loading Rack (Set of 2)	Set of 2 racks, each accommodates 23 reagent vials. Adapters not included	662664
Reagent Rack Adapters		
Reagent Rack Adapter Type 1 (Pack of 10)	Adapter to accommodate BD TruCount™ controls, and 4mL BD Horizon™ Brilliant Stain Buffer	663260
Reagent Rack Adapter Type 2 (Pack of 10)	Adapter to accommodate 0.5mL and 2mL Pharmingen Reagent Vials	663261
Miscellaneous Accessories		
Lysing Solution 1L Tank	1L tank for lysing solution	662602
10L DI Water Tank	10L tank for DI water	662603
10L Waste Tank	10L tank for waste	662604
10L Saline Tank	10L tank for saline	662605
50mL Bleach Container	50mL bleach container for end of day clean cycle	662606
50mL Saline Container	50mL container for saline (in the reagent bay)	662607

Maintenance kits

Item	Description	Cat. No.
Peristaltic Pump Tubing Replacement Kit	Kit for replacement of the tubing for peristaltic pumps	662613
Specimen Probe Replacement Kit	Kit for replacement of specimen pipettor probe	662614
Syringe Pump Valve Replacement Kit	Kit for replacement of the syringe valve pump	663247
Syringe Pump Replacement Kit	Kit for replacement of the 500µL syringe	663249
DI Water Filter Replacement Kit	Kit for replacement of the DI water filter	663923
Maintenance kit	Kit for 6 month maintenance (includes 662614, 662613, 663247, 663249 and 663923)	663926

System and LIS connectivity

Item	Description	Cat. No.
BD FACSLink™ Basic Software + 1 Year Maintenance (EU only)	BD FACSLink basic software + 2 connections + 1 year maintenance	663061
BD FACSLink Computer Assembly (EU only)	Computer assembly for BD FACSLink software	663874
BD FACSLink User License + 1 Year Maintenance (EU only)	Extra connection license (optional, for multiple instruments)	663062
BD FACSLink Computer Assembly (non-EU)	Computer assembly for BD FACSLink software	663875
BD FACSLink Software Package (non-EU)	BD FACSLink software for BD FACSLyric. Allows for connection to LIS	663893
BD FACSLink Connection License for Additional Instruments (non-EU)	Extra BD FACSLink connection license (optional, for multiple instruments)	663625
BD FACSLink Software Media Kit (for all regions except EU, US, Canada)	Media kit required for first time installations or software upgrades of BD FACSLink	663727

BD FACSDuet[™] is for In Vitro Diagnostic Use. CE marked in compliance with the European In Vitro Diagnostic Medical Device Directive 98/79/EC.

X23-20471-01EU

BD Life Sciences, San Jose, CA, 95131, USA

bdbiosciences.com/eu