## **Maintenance Log for Sample Preparation Instruments**

																			ľ	<b>llo</b> n	th/	<b>Yea</b>	r: _								
Serial No.:	Initials and Dates																														
SAMPLE PREP																															
ASSISTANT III MAINT	ENA	NCE																													
	Da	ily																													
DATE																										Ī	Ī				
Daily Clean																															
	We	ekly														,					•										
Weekly Clean																															
Calibrate Probe Reference Point																															
	Mo	nthly	•							·															·						
Waste In-Line Filter Service																															
	Pro	be C	hang	је	I					<u> </u>					1																
Change Probe																															
Verify Accuracy and Precision																															

Serial No.:	Initials and Dates																					
LYSE WASH ASSISTA	NT																					
MAINTENANCE																						
	Daily	/																				
DATE																						
Daily Wash																						
	Week	dy	•							•	•				.,				•	1		•
Clean the Cell Wash Station																						
	Mont	hly												ı								
Change waste tank cap																						
	Bi-An	nual																				
Clean the Fluidics Tanks																						
Replace the Bal Seal																						
	ı								1					1								