

BD Accuri™ C6 Plus System

Product list



BD Accuri[™] C6 Plus System

Flow cytometry systems

Description	Notes	Cat. no.
BD Accuri™ C6 Plus Flow Cytometer System	Includes standard optical filters, reagent containers, complimentary BD Accuri™ C6 Plus Flow Cytometer starter kit with 2-month maintenance kit and 1-year warranty	660517
BD Accuri™ C6 Plus Workstation Computer and Software	Computer and software required to acquire and analyze samples from the BD Accuri™ C6 Plus Flow Cytometer	663931
BD CSampler™ Plus: 24-Tube Rack and 96-Well Plate Automated Sample Loader	Includes 24-tube rack (white) and 1-year warranty for BD Accuri™ C6 Plus Flow Cytometer	660519

Optional system modules

Description	Notes	Cat. no.
Selectable lasers module for BD Accuri™ C6 Plus Flow Cytometer	Includes activation key and optical filters: 780/60, 610/20, 630/30 Facilitates standard three detectors off blue laser, and one off red laser Allows for four blue detectors, or two blue and two red detectors	661085
User tracking module for BD Accuri™ C6 Plus Flow Cytometer	Includes activation key. Enables the creation of unique user names and passwords to record the duration that each operator uses the cytometer	661086
Remote interfacing module for BD Accuri™ C6 Plus Flow Cytometer	Includes activation keys. Enables BD Accuri™ C6 Plus or BD CSampler™ Plus Software to receive commands from other software applications via a TCP/IP connection	661087

Flow cytometry analysis software

Description	Notes	Cat. no.
BD Accuri™ C6 Plus Analysis Software for PC or Mac	For use with data exported from the BD Accuri™ C6 Plus Flow Cytometer	664022
BD CSampler™ Plus Analysis Software for PC or Mac	For use with data exported from the BD Accuri™ C6 Plus Flow Cytometer with BD CSampler™ Plus Software	664023
FCS Express™ 5 Flow Cytometry Software, RUO, single license	With BD Accuri™ C6 Plus File Importer	645024

Extended warranty (purchased after instrument purchase)

Description	Notes	Cat. no.
Service contract for BD Accuri™ C6 Plus Flow Cytometer (per year)	For purchasing 1-year or multiple-year service contract on BD Accuri™ C6 Plus Flow Cytometer	661498
Service contract for BD Accuri™ Plus Flow Cytometer (per year)	For purchasing 1-year or multiple-year service contract on BD CSampler™ Plus Rack	661499

Maintenance, quality assurance and service

Description	Notes	Cat. no.
BD® CS&T RUO Beads (50 tests)	For use with BD Accuri™ C6 Plus and BD Accuri™ Plus Software	661414
BD® CS&T RUO Beads (150 tests)	For use with BD Accuri™ C6 Plus and BD Accuri™ Plus Software	661415
Peristaltic pump tubing (set of 2)		659522
Sheath, detergent, BD FACSClean™ Solution Bottle Filters (1 each)		660322
In-line sheath filter (1)		660499
BD Accuri™ C6 Plus Flow Cytometer 1-year maintenance kit	Contains peristaltic pump tubing (12), sheath bottle filters (6) and in-line sheath filters (6)	661331
BD Accuri™ C6 Plus Flow Cytometer 2-month maintenance kit	Contains peristaltic pump tubing (2), sheath bottle filters (1) and in-line sheath filters (1)	660634
BD Accuri™ C6 Plus Flow Cytometer install and training		661269
BD Accuri™ C6 Plus Flow Cytometer onsite training		661290
BD Accuri™ C6 Plus Flow Cytometer IQ/OQ		661526

BD Accuri™ C6 Plus System

Reagents and kits

Description	Notes	Cat. no.
BD Accuri™ C6 Plus Flow Cytometer Starter Kit	Contains concentrate for making 1 liter each of BD® Detergent Solution and BD FACSClean™ Solution, BD® Sheath Additive concentrate for treating 10 L of sheath, 8 mL of our BD® Extended Flow Cell Clean Solution, BD® CS&T RUO Beads for 50 tests, peristaltic pump tubing (set of 2)	661535
BD Accuri™ C6 Plus Flow Cytometer Fluidic Kit	Contains concentrate for making 1 liter each of BD® Detergent Solution and ter Fluidic Kit BD FACSClean™ Solution, BD® Sheath Additive concentrate for treating 10 L of sheath, 8 mL of our BD® Extended Flow Cell Clean Solution	
BD® Sheath Additive		660584
BD® Detergent Solution Concentrate		660585
BD FACSClean™ Solution		340345
BD® Extended Flow Cell Clean Solution		660586

Replacement parts and accessories

Description	Notes	Cat. no.
BD Accuri™ Flow Cytometer Protector Pelican Case	Includes packaging foam	653163
BD CSampler™ Protector Pelican Case	Includes packaging foam	653162
24-tube rack with key		661274
BD Accuri™ 24-Tube Rack		653161
2 L waste bottle		659605
2 L sheαth bottle		661277
Detergent bottle		661276
BD FACSClean™ Bottle		661275

Optical filters

Description	Cat. no.
FL1 99% attenuation filter	653172
FL1 90% attenuation filter	653173
FL2 99% attenuation filter	653174
FL2 90% attenuation filter	653175
FL3 99% attenuation filter	653176
FL3 90% attenuation filter	653177
FL4 99% attenuation filter	653178
FL4 90% attenuation filter	653179
Standard optical filter 533/30 nm (FL1)	653180
Standard optical filter 585/40 nm (FL2)	653181
Standard optical filter >670 nm LP (FL3)	653182
Standard optical filter 675/25 nm (FL4)	653183
Optical filter 510/15 nm	653184
Optical filter 565/20 nm	653185
Optical filter 610/20 nm	653186
Optical filter 780/60 nm	653187
Optical filter 540/20 nm	653528

Please note flow cytometer serial number prior to contacting customer service to ensure purchase of correct optical filters.

Learn more about our products at **bdbiosciences.com**

Applicable terms and conditions are available at bdbiosciences.com.
Class 1 Laser Product.
For Research Use Only. Not for use in diagnostic or therapeutic procedures.
FCS Express is a trademark of De Novo Software.
23-18355-02

BD, San Jose, CA, 95131, U.S.

bd.com

