



BD FACSDuet™ Sample Preparation System

Product list



BD FACSDuet™ Sample Preparation System

Sample Preparation System

Item	Description	Cat. No.
BD FACSDuet™ Sample Preparation System	Sample preparation system for all Lyse-No-Wash applications. Includes kit for physical integration with the BD FACSLytic™ flow cytometer. Also includes a starter accessory kit:	662588
	Reagent Loading Rack (Set of 2 racks)	
	Specimen Loading Rack (Set of 4 racks)	
	13x75 Specimen Tube Adapters (4 packs of 10 ea.)	
	BD FACSLytic Secondary Tube Rack (30 Tube Rack) (Set of 4 racks)	
	Handheld Barcode Scanner with Stand	
	Lysing Solution 1L Tank	
	10L DI Water Tank	
	10L Waste Tank	
	10L Saline Tank	
	50mL Bleach Container	
	50mL Chilled Saline Container	
	Documentation Kit (IFU, User Guide, SPG, SLG)	
	Panel PC	

Optional components

Item	Description	Cat. No.
Table (recommended when BD FACSDuet and BD FACSLytic are a physically integrated system)	200.0cm L x 97.0cm D x 86.35cm H Table includes BD FACSLytic flat screen desk support, dripping plate with receptacles, 2 chariots for BD FACSDuet and BD FACSLytic tanks	663375
Secondary Tube Labels Printer	Printer generating labels for 12x75mm secondary tubes	662609
Paper for secondary tube label printer	3 rolls + 1 ribbon	659170
2D Barcode Zebra printer Model ZD410	Printer generating 2D barcode for reagent vial labels	257705
Labels roll for reagent vial labels	1500 +/- 50 labels; label size 2" x 0.75"	441404
BD Unlabeled Empty Reagent Amber Vials	1 box with 24 glass vials + cap	653782
BD 2D Barcoded Empty Reagent Amber Vials	1 box with 24 glass vials + cap	664358
Handheld Barcode Scanner, with Stand	External barcode reader for both BD FACSLytic and BD FACSDuet	663058
Uninterruptible Power Manager (UPM)	1440VA 230V IN/OUT	663407

BD FACSDuet™ Sample Preparation System

Accessories


Item	Description	Cat. No.
Specimen Loading Rack		
Specimen Loading Rack (Set of 4)	Set of 4 racks for loading primary specimen tubes. Does not include tube adaptors	662592
Streck Cyto-Chex specimen Loading Rack (Set of 4)	Set of 4 racks for loading Streck Cyto-Chex tubes.	662593
Specimen Rack Adapters		
8x66 Specimen Tube Adapters (set of 10)	8mmx66mm adapter for Sarstedt 1.2 mL tubes	664848
13x65 Specimen Tube Adapters (Set of 10)	13mmx65mm adapter for Sarstedt 2.6mL and 3.4mL tubes	662594
11x66 Specimen Tube Adapters (Set of 10)	11mmx66mm adapter for Sarstedt 2.7mL tubes	662595
13x75 Specimen Tube Adapters (Set of 10)	13mmx75mm adapter for BD 2mL, 3mL and 4mL, Greiner 2mL and 3mL, Streck 5mL tubes	662596
15x75 Specimen Tube Adapters (Set of 10)	15mmx75mm adapter for Sarstedt 4mL tubes	662597
13x90 Specimen Tube Adapters (Set of 10)	13mmx90mm adapter for Sarstedt 4.9mL tubes	662598
13x100 Specimen Tube Adapters (Set of 10)	13mmx100mm adapter for BD and Greiner 6mL and 7mL tubes	662600
16x100 Specimen Tube Adapters (Set of 10)	16mmx100mm adapter for BD 9.5mL and 10mL and Greiner 9mL and 10mL tubes	662599
Secondary Tube Racks		
BD FACSLyric Secondary Tube Rack (30 Tube Rack)	Secondary tube rack to accommodate 30 12x75mm tubes	651314
BD FACSLyric Secondary Tube Rack (40 Tube Rack)	Secondary tube rack to accommodate 40 12x75mm tubes	651319
Reagent Racks		
Reagent Loading Rack (Set of 2)	Set of 2 racks, each accommodates 23 reagent vials. Adapters not included	662664
Dual-size Reagent Loading Rack (Set of 1)	Set of 1 rack, accommodates 23 reagent vials - 12 Beckman Coulter, 11 for BD/others. Adapters not included.	664458
Reagent Vial Adapters		
Reagent Vial Adapter Type 1 (Pack of 10)	Adapter to accommodate BD TruCount™ controls, and 4mL BD Horizon™ Brilliant Stain Buffer	663260
Reagent Vial Adapter Type 2 (Pack of 10)	Adapter to accommodate 0.5mL and 2mL Pharmingen Reagent Vials	663261
Reagent Vial Adapter Type 3 (Pack of 5)	Adapter to accommodate Cytognos 4-mL reagent vials	664457
Reagent Vial Adapter Type 4 (Pack of 5)	Adapter to accommodate Beckman Coulter 1-mL reagent vials - requires Dual-size Reagent Loading Rack (664458)	664459
Reagent Vial Adapter Type 5 (Pack of 5)	Adapter to accommodate Dako 6-mL reagent vials	664460
Miscellaneous Accessories		
Lysing Solution 1L Tank	1L tank for lysing solution	662602
10L DI Water Tank	10L tank for DI water	662603
10L Waste Tank	10L tank for waste	662604
10L Saline Tank	10L tank for saline	662605
50mL Bleach Container	50mL bleach container for end of day clean cycle	662606
50mL Saline Container	50mL container for saline (in the reagent bay)	662607

Maintenance kits

Item	Description	Cat. No.
Peristaltic Pump Tubing Replacement Kit	Kit for replacement of the tubing for peristaltic pumps	662613
Specimen Probe Replacement Kit	Kit for replacement of specimen pipettor probe	662614
Syringe Pump Valve Replacement Kit	Kit for replacement of the syringe valve pump	663247
Syringe Pump Replacement Kit	Kit for replacement of the 500µL syringe	663249
DI Water Filter Replacement Kit	Kit for replacement of the DI water filter	663923
Maintenance kit	Kit for 6 month maintenance (includes 662614, 662613, 663247, 663249 and 663923)	663926

System and LIS connectivity

Item	Description	Cat. No.
BD FACSLink™ Basic Software + 1 Year Maintenance (EU only)	BD FACSLink basic software + 2 connections + 1 year maintenance	663061
BD FACSLink Computer Assembly (EU only)	Computer assembly for BD FACSLink software	664487
BD FACSLink User License + 1 Year Maintenance (EU only)	Extra connection license (optional, for multiple instruments)	663062
BD FACSLink Computer Assembly (non-EU)	Computer assembly for BD FACSLink software	663875
BD FACSLink Software Package (non-EU)	BD FACSLink software for BD FACSLyric. Allows for connection to LIS	663893
BD FACSLink Connection License for Additional Instruments (non-EU)	Extra BD FACSLink connection license (optional, for multiple instruments)	663625
BD FACSLink Software Media Kit (for all regions except EU, US, Canada)	Media kit required for first time installations or software upgrades of BD FACSLink	663727




BD FACSDuet™ Sample Preparation System and the BD FACSLyric™ Flow Cytometer with the BD FACSuite™ Clinical and BD FACSuite™ applications are CE marked in compliance with the European In Vitro Diagnostic Medical Device Directive 98/79/EC.

BD FACSDuet™ Sample Preparation System and BD FACSLyric™ Flow Cytometer are Class 1 Laser Products.

BD Switzerland Sàrl, Terre-Bonne Park - A4, Route de Crassier -17, 1262 Eysins, Switzerland

bdbiosciences.com/eu



 BD, the BD Logo, FACSDuet, FACSLink, FACSLyric and FACSuite are trademarks of Becton, Dickinson and Company or its affiliates. © 2021 BD. All rights reserved. X23-20471-03.