Technical Data Sheet

Purified Mouse Anti-Human CD209

Product Information

Material Number:551186Alternate Name:DC-SIGNSize:0.1 mgConcentration:0.5 mg/mlClone:DCN46


Immunogen: Human Monocyte Derived DC Cells

Isotype:Mouse IgG2b, κ Reactivity:QC Testing: Human

Storage Buffer: Aqueous buffered solution containing ≤0.09% sodium azide.

Description

The DCN46 antibody specifically binds to dendritic cell-specific ICAM-3 grabbing nonintegrin (DC-SIGN or CD209), a type-II membrane protein of approximately 44 kDa with a mannose-binding C-type lectin domain. It is highly expressed on dendritic cells in mucosal tissues. Its sequence is identical to the HIV-1 envelope gp120-binding C-type lectin, and reports suggest that DC-SIGN binds to HIV-1 gp120 and effectively transmits infectious HIV-1 to resting T lymphocytes expressing CD4 and chemokine receptors. The C-type lectin domain of DC-SIGN is also capable of binding other pathogenic viruses, bacteria, and parasites. Reports also suggest that DC-SIGN enables the highly efficient migration of dendritic cells from blood into the tissues. It can interact with ICAM-2, which has a similar sequence as ICAM-3, and is abundantly expressed on vascular and lymphoid endothelium. Thus, DC-SIGN mediates dendritic cells rolling and transendothelial migration, and its interaction with ICAM-2 is essential to specific migratory functions of dendritic cells.



Flow cytometric analysis of CD209 expression on cultured human monocyte-derived dendritic cells.
Dendritic cells were stained with either Purified Mouse
Anti-Human CD209 (Cat. No. 551186; solid line histogram)
or Purified Mouse IgG2b κ Isotype Control (Cat. No. 555740;
dashed line histogram), followed by FITC Goat Anti-Mouse
IgG/IgM (Cat. No. 555988). Fluorescent histograms were
derived from gated events with the side and forward
light-scattering characteristics of viable cells.

Preparation and Storage

Store undiluted at 4°C.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

Application Notes

Application

| /1000VI | | |
|----------------|------------------|--|
| Flow cytometry | Routinely Tested | |

Suggested Companion Products

| Catalog Number | Name | Size | Clone |
|----------------|--|--------|------------|
| 555740 | Purified Mouse IgG2b κ Isotype Control | 0.1 mg | 27-35 |
| 555988 | FITC Goat Anti-Mouse IgG/IgM | 0.5 mg | Polyclonal |
| 554657 | Stain Buffer (BSA) | 500 mL | (none) |
| 554656 | Stain Buffer (FBS) | 500 mL | (none) |

BD Biosciences

bdbiosciences.com

United States Canada Europe Japan Asia Pacific Latin America/Caribbean 877.232.8995 866.979.9408 32.2.400.98.95 0120.8555.90 65.6861.0633 55.11.5185.9995

For country contact information, visit ${\bf bdbiosciences.com/contact}$

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton, Dickinson and Company is sticity prohibited.

For Research Use Only, Not for use in diagnostic or therapeutic procedures. Not for resale.

© 2017 BD. BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.



551186 Rev. 9 Page 1 of 2

Product Notices

- 1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
- 2. An isotype control should be used at the same concentration as the antibody of interest.
- 3. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
- 4. Sodium azide is a reversible inhibitor of oxidative metabolism; therefore, antibody preparations containing this preservative agent must not be used in cell cultures nor injected into animals. Sodium azide may be removed by washing stained cells or plate-bound antibody or dialyzing soluble antibody in sodium azide-free buffer. Since endotoxin may also affect the results of functional studies, we recommend the NA/LE (No Azide/Low Endotoxin) antibody format, if available, for in vitro and in vivo use.
- 5. Please refer to www.bdbiosciences.com/pharmingen/protocols for technical protocols.

References

Geijtenbeek TBH, Kwon DS, Torensma R, et al. DC-SIGN, a dendritic cell-specific HIV-1 -binding protein that enhances trans-infection of T cells. Cell. 2000; 100(5):587-597. (Biology)

Geijtenbeek TBH, Torensma R, van Vliet SJ, et al. Identification of DC-SIGN, a novel dendritic cell-specific ICAM-3 receptor that supports primary immune responses. Cell. 2000; 100(5):575-585. (Biology)

Sallusto F, Cella M, Danieli C, Lanzavecchia A. Dendritic cells use macropinocytosis and the mannose receptor to concentrate macromolecules in the major histocompatibility complex class II compartment: downregulation by cytokines and bacterial products. *J Exp Med.* 1995; 182(2):389-400. (Immunogen) Steinman RM. DC-SIGN: a guide to some mysteries to dendritic cells. *Cell.* 2000; 100(5):491-494. (Biology)

551186 Rev. 9 Page 2 of 2