

4.5 Ounce Container, with and without lid - Total Volume = 110 ml

1.0 Applications:

Collection , transportation, and storage of a wide variety of specimens

2.0 Features:

- Sterile, disposable containers.
- Molded of rugged polypropylene
- Inert and chemically resistant to commonly used laboratory reagents at room temperature.
- Provide secure, time saving sample containment.
- Molded-in graduations for easy measurements.
 - Graduated in ounces and milliliters.
 - Graduated from 1/2 oz. to 4.5 oz. increments and from 20mL to 110mL in 10 mL increments.
- Supplied with or without polyethylene snap lid.

3.0 STERILITY:

Product is gamma irradiated and dosimetrically released based upon U.S. Association for the Advancement of Medical Instrumentation (AAMI) recommended practices in effect at the time of validation. Sterilization records are reviewed and signed off by qualified personnel for product release.

4.0 FDA COMPLIANCE:

FALCON products are manufactured in accordance with the current FDA Quality System Regulation (1996) as specified in Chapter 1, Title 21 of the Code of USA Federal Regulations, Part 820 for medical devices.

4.0 RESTRICTIONS:

Refer to "Thermoplastics Properties Chart," available via internet at <http://www.bd.com/labware/freestuf/>

4.5 Ounce Container, with and without lid - Total Volume = 110 ml

Related Catalog Numbers:


Cat. No.	Packout Description	Lid	Sterile	Graduations
354013	1/Bag, 100/Case	Yes	Yes	Molded
354014	20/Bag, 500/Case	No	Yes	Molded

Dimensions: XXXmm (XXXinches)

LID: Catalog No. 354017 20/Bag, 500/Case

Material: Polyethylene

A = O.D. 68.81 (2.709)
B = O.D. 71.12 (2.800)
C = 7.75 (.305)



CONTAINER:

Material: Polypropylene

A = O.D. 67.18 (2.645)
B = O.D. 52.30 (2.059)
C = 64.21 (2.528)


Top I.D. below rim = 61.90 (2.437)

Bottom I.D. = 51.49 (2.027)

Wall Thickness

Side Wall = .48 (.019)

Bottom Wall = .56 (.022)



Container Height w/Lid = 65.41 (2.575)