

Technical Data Sheet


Purified NA/LE Mouse Anti-Human CD28

Product Information

Material Number:	555725
Alternate Name:	CD28 antigen; T44; Tp44; TP44
Size:	0.5 mg
Concentration:	1.0 mg/ml
Clone:	CD28.2
Isotype:	Mouse IgG1, κ
Reactivity:	QC Testing: Human
Workshop:	V 5T CD28.05
Storage Buffer:	No azide/low endotoxin: Aqueous buffered solution containing no preservative, 0.2 μ m sterile filtered. Endotoxin level is \leq 0.1 EU/ μ g (\leq 0.001 ng/ μ g) of protein as determined by the LAL assay.

Description

The CD28.2 monoclonal antibody specifically binds to CD28, a 44 kDa homodimeric transmembrane glycoprotein present on most mature T cells, thymocytes and plasma cells. CD28 is a costimulatory receptor that binds CD80 and CD86 as ligands and plays a very important role in T cell-B cell interactions. It has been suggested that CD28 initiates and regulates a separate and distinct signal transduction pathway from those stimulated by the TCR complex. Additionally, it has been reported that CD28 antibody clones vary in their ability to stimulate T cells to produce IL-2 and increase intracellular Ca²⁺ concentration. This finding suggests the existence of functionally distinct subregions on the CD28 molecule. CD28.2 has been demonstrated to bind to the same molecule as clone L293, another CD28 mAb, and has been reported to induce Ca²⁺ influx in Jurkat T cells.



Profile of peripheral blood lymphocytes analyzed by flow cytometry. Second step staining with Cat. No. 555988.

Preparation and Storage

Store undiluted at 4°C.

This preparation contains no preservatives, thus it should be handled under aseptic conditions.

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography.

Application Notes

Application

Flow cytometry	Routinely Tested
Bioassay	Tested During Development

Recommended Assay Procedure:

Suggested protocol for functional studies using CD28.2 mAb:

1. Isolate human PBMC with Ficoll-Paque™ Plus
2. Suspend PBMC in medium supplemented with 1% glutamine, 1% penicillin/streptomycin and 10% FBS at 10⁶ cells/ml
3. Incubate 1 x 10⁶ PBMC with soluble NA/LE™ format of CD3 mAb (Cat. No. 555336) at 1 μ g/ml final concentration in 7% CO₂, 37°C for 3

BD Biosciences

bdbiosciences.com

United States	Canada	Europe	Japan	Asia Pacific	Latin America/Caribbean
877.232.8995	888.259.0187	32.53.720.550	0120.8555.90	65.6861.0633	55.11.5185.9995

For country-specific contact information, visit bdbiosciences.com/how_to_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2008 BD



days

4. Wash cells twice and resuspend in the above medium at 10^6 cells/ml
5. Distribute cells in 96-well round bottom plate at 10^5 cells/well
6. Add soluble NA/LE™ format of CD28.2 mAb at 5 µg/ml final concentration
7. Incubate the plate for an additional three days
8. Pulse the plate with [³H]TdR (1 µCi/well). Harvest the plate 20 hours after pulse and count

Comment: Addition of Protein G (Sigma) 5 µg/ml together with CD28.2 will significantly enhance the costimulation, presumably through cross-linking.

Suggested Companion Products

<u>Catalog Number</u>	<u>Name</u>	<u>Size</u>	<u>Clone</u>
555988	FITC Goat Anti-Mouse IgG/IgM	0.5 mg	Polyclonal
555336	Purified NA/LE Mouse Anti-Human CD3	0.5 mg	HIT3a

Product Notices

1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
2. Ficoll-Paque is a trademark of Amersham Biosciences Limited.
3. Please refer to www.bdbiosciences.com/pharming/en/protocols for technical protocols.

References

- Barclay NA, Brown MH, Birkeland ML, et al, ed. *The Leukocyte Antigen FactsBook*. San Diego, CA: Academic Press; 1997. (Biology)
- Chen Z. Personal Communication. . (Biology)
- June CH, Bluestone JA, Nadler LM, Thompson CB. The B7 and CD28 receptor families. *Immunol Today*. 1994; 15(7):321-331. (Biology)
- Nunes J, Klasen S, Franco MD, et al. Signalling through CD28 T-cell activation pathway involves an inositol phospholipid-specific phospholipase C activity. *Biochem J*. 1993; 293(3):835-842. (Biology)
- Nunes J, Klasen S, Ragueneau M, et al. CD28 mAbs with distinct binding properties differ in their ability to induce T cell activation: analysis of early and late activation events. *Int Immunol*. 1993; 5(3):311-315. (Biology)
- Schlossman SF, Boumsell L, Gilks W, et al, ed. *Leukocyte Typing V: White Cell Differentiation Antigens*. New York: Oxford University Press; 1995. (Clone-specific)
- Verwilghen J, Vandenberghe P, Wallays G, et al. Simultaneous ligation of CD5 and CD28 on resting T lymphocytes induces T cell activation in the absence of T cell receptor/CD3 occupancy. *J Immunol*. 1993; 150(3):835-846. (Biology)