

## Technical Data Sheet

**Purified Mouse Anti-Fyn (pY528)/c-Src (pY530)****Product Information**

<b>Material Number:</b>	612669
<b>Size:</b>	150 µg
<b>Concentration:</b>	250 µg/ml
<b>Clone:</b>	31/Fyn (pY528)/c-Src (pY530)
<b>Immunogen:</b>	Phosphorylated Human Fyn (Y528) Peptide
<b>Isotype:</b>	Mouse IgG2b
<b>Reactivity:</b>	QC Testing: Human Tested in Development: Mouse
<b>Target MW:</b>	60 kDa
<b>Storage Buffer:</b>	Aqueous buffered solution containing BSA, glycerol, and ≤0.09% sodium azide.

**Description**

p59fyn is a member of the Src family of protein tyrosine kinases. Two isoforms of the protein have been identified, designated as FynB and Fyn T. Fyn B has been shown to be localized to the brain, whereas Fyn T associates with both B and T cells. Stimulation of the T cell antigen receptor (TcR) results in protein tyrosine phosphorylation via non-receptor tyrosine kinases. p59fyn kinase (Fyn T) associates with the TcR. Ligation of the TcR activates the protein kinase activity of p59fyn in various human T cells. Fyn interacts with the CD3-zeta chains through its N-terminal region. In turn, fyn binds other proteins through its SH2 and SH3 domains. These proteins (p82 and p116) may serve as substrates and/or mediators of fyn activity. p59fyn kinase is tyrosine phosphorylated at two sites, Tyr-417 (autophosphorylation site), and Tyr-528 (negative regulatory site).



**Western blot analysis for Fyn (pY528) / c-Src (pY530).**  
A431 cells (Human epithelial carcinoma; ATCC CRL-1555) were treated with 100 ng/ml EGF for 5 min and then either left untreated (lane 1) or treated (lane 2) with 50 µg/ml alkaline phosphatase for 30 min at 37°C. The top panel was probed with a mouse anti-Src antibody (Millipore/Upstate Cat. No. 05-184) and the bottom was probed with the mouse anti-Fyn (pY528) / c-Src (pY530) antibody at a 1:2500 dilution with a band observable at ~ 60 kDa.

**Preparation and Storage**

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography. Store undiluted at -20°C.

**Application Notes****Application**

Western blot	Routinely Tested
--------------	------------------

**Recommended Assay Procedure:**

**Western blot:** Please refer to [http://www.bdbiosciences.com/pharmingen/protocols/Western\\_Blotting.shtml](http://www.bdbiosciences.com/pharmingen/protocols/Western_Blotting.shtml)

**Suggested Companion Products**

Catalog Number	Name	Size	Clone
611448	A431 + EGF Cell Lysate	500 µg	(none)
554002	HRP Goat Anti-Mouse Ig	1.0 ml	(none)

**BD Biosciences**

bdbiosciences.com

United States	Canada	Europe	Japan	Asia Pacific	Latin America/Caribbean
877.232.8995	888.259.0187	32.53.720.550	0120.8555.90	65.6861.0633	55.11.5185.9995

For country-specific contact information, visit [bdbiosciences.com/how\\_to\\_order/](http://bdbiosciences.com/how_to_order/)

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2008 BD



## Product Notices

1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
2. Please refer to [www.bdbiosciences.com/pharmingen/protocols](http://www.bdbiosciences.com/pharmingen/protocols) for technical protocols.
3. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.
4. Source of all serum proteins is from USDA inspected abattoirs located in the United States.

## References

- Cooke MP, Perlmutter RM. Expression of a novel form of the fyn proto-oncogene in hematopoietic cells. *New Biol.* 1989; 1(1):66-74.(Biology)
- Grant SG, O'Dell TJ, Karl KA, Stein PL, Soriano P, Kandel ER. Impaired long-term potentiation, spatial learning, and hippocampal development in fyn mutant mice. *Science.* 1992; 258(5090):1903-1910.(Biology)
- Kawakami T, Kawakami Y, Aaronson SA, Robbins KC. Acquisition of transforming properties by FYN, a normal SRC-related human gene. 1988; 85(11):3870-3874.(Biology)