

Technical Data Sheet

Purified Mouse Anti-Rat CD45RA

Product Information

| | |
|-------------------------|--|
| Material Number: | 554882 |
| Size: | 0.5 mg |
| Concentration: | 0.5 mg/ml |
| Clone: | OX-33 |
| Immunogen: | Leukocyte common antigen purified from rat splenocytes |
| Isotype: | Mouse (BALB/c) IgG1, κ |
| Reactivity: | QC Testing: Rat |
| Storage Buffer: | Aqueous buffered solution containing $\leq 0.09\%$ sodium azide. |

Description

The OX-33 antibody reacts with a high-molecular-weight form of CD45 found only on B lymphocytes. CD45 is a member of the Protein Tyrosine Phosphatase (PTP) family. Its intracellular (COOH-terminal) region contains two PTP catalytic domains, and the extracellular region is highly variable due to alternative splicing of exons 4, 5, and 6 (designated A, B, and C, respectively), plus differing levels of glycosylation. The CD45 isoforms detected in the rat are cell type-, maturation-, and activation state-specific. The CD45 isoforms play complex roles in T-cell and B-cell antigen receptor signal transduction.

This antibody is routinely tested by flow cytometric analysis. Other applications were tested at BD Biosciences Pharmingen during antibody development only or reported in the literature.

Preparation and Storage

The monoclonal antibody was purified from tissue culture supernatant or ascites by affinity chromatography. Store undiluted at 4° C.

Application Notes

Application

| | |
|---------------------------------|------------------|
| Flow cytometry | Routinely Tested |
| Immunohistochemistry-zinc-fixed | Reported |
| Immunohistochemistry-frozen | Reported |

Recommended Assay Procedure:

We have observed that the staining intensity of OX-33 mAb is reduced after fixation of stained leukocytes (≤ 3 hours with 1% formaldehyde). Therefore, one should not fix the stained cells prior to flow cytometry. We have found that freshly-isolated leukocytes and cell lines may wait for analysis in wash buffer at 4°C, without fixation, for up to 18 hours post-staining without loss of viability. Activated lymphocytes may lose viability rapidly, and data should be collected within 5 hours post-staining. Other reported applications include immunohistochemical staining (IHC) of acetone-fixed frozen sections. mAb OX-33 is not recommended for IHC of formalin-fixed paraffin-embedded sections.

Suggested Companion Products

| Catalog Number | Name | Size | Clone |
|----------------|---|--------|------------|
| 555988 | FITC Goat Anti-Mouse IgG/IgM | 0.5 mg | Polyclonal |
| 557273 | Purified Mouse IgG1, κ Isotype Control | 0.5 mg | MOPC-31C |

Product Notices

1. Since applications vary, each investigator should titrate the reagent to obtain optimal results.
2. Please refer to wwwbdbiosciences.com/pharmingen/protocols for technical protocols.
3. Caution: Sodium azide yields highly toxic hydrazoic acid under acidic conditions. Dilute azide compounds in running water before discarding to avoid accumulation of potentially explosive deposits in plumbing.

BD Biosciences

bdbiosciences.com

United States 877.232.8995 Canada 888.259.0187 Europe 32.53.720.550 Japan 0120.8555.90 Asia Pacific 65.6861.0633 Latin America/Caribbean 55.11.5185.9995

For country-specific contact information, visit bdbiosciences.com/how_to_order/

Conditions: The information disclosed herein is not to be construed as a recommendation to use the above product in violation of any patents. BD Biosciences will not be held responsible for patent infringement or other violations that may occur with the use of our products. Purchase does not include or carry any right to resell or transfer this product either as a stand-alone product or as a component of another product. Any use of this product other than the permitted use without the express written authorization of Becton Dickinson and Company is strictly prohibited.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Not for resale.

BD, BD Logo and all other trademarks are the property of Becton, Dickinson and Company. ©2006 BD



BD

BD Biosciences

References

- Johnson P, Maiti A. CD45: A family of leukocyte-specific cell surface glycoproteins. In: Herzenberg LA, Weir DM, Blackwell C, ed. *Weir's Handbook of Experimental Immunology, Vol 2*. Cambridge: Blackwell Science; 1997:62.1-62.16.(Biology)
- Woollett GR, Barclay AN, et al. Molecular and antigenic heterogeneity of the rat leukocyte common antigen from thymocytes and T and B lymphocytes.. *Eur J Immunol*. 1985; 15:168-173.(Immunogen: Flow cytometry)